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DRY INCUBATOR- USER MANUAL

Remarks:

- ◆ These operator's instructions have to be carefully kept near the equipment, in order to be consulted by all users.
- ◆ In case of lack of voltage or lighting of the alarm pilot lamp the equipment stops the incubation cycle.
- ◆ We recommend not to lean any object on the equipment.
- ◆ Temperature oscillations of few °C, during the cycle, have to be considered normal.
- ◆ Always place one of the pierced grids on the lower plane (1).

Description of the equipment :

Referring to the figure of page 2 the incubator components are described here below:

- **(A) Master switch**
- **(B) "Ready machine" yellow pilot lamp**

This lamp is on when the incubator is ready and is waiting the start of a cycle

- **(C) Alarm red lamp** 

This pilot lamp is turned on in order to signal some possible anomalies (for ex. lack of voltage). When the master switch is on, the pilot lamp is turned on to signal the previous switching off. By pressing the START key the alarm is cancelled. If, during the cycle, this pilot lamp is turned on, the cycle is to be considered null (see also Alarms Table).

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The function of this push button is to start a cycle, if the pilot lamp (B) is turned on.

The **START** push button is also working as RESET key, if it is pressed during a incubation cycle the cycle is cancelled. Moreover, if an alarm is intervening, by pressing the **START** push button, the alarm is cancelled.

- **(E) Switch push button** 

The function of this push button is to alternate, on the display (H), the visualization of temperature and incubation time.

- **(F) Temperature pilot lamp** 

When this pilot lamp is turned on it shows that the value read and/or set on the display (H) is referred to the temperature.

- **(G) time pilot lamp** $\frac{1}{2}$

When this pilot lamp is turned on it shows that the value read and/or set on the display (H) is referred to the incubation length of time, heating excluded.

- **(H) display**

Cycle parameters are visualized on the display : temperature inside the chamber, incubation time, set temperature.

In order to reset the cycle data, this must be stopped.

- The temperature can be set between 30°C and 70°C,

If the temperature inside the chamber is lower than 25°C you can read "t.Lo", (low temperature) on the display, if the temperature is higher than 110°C you can read "t.Hi" (high temperature).

- Time can be set between 1 minute and 120 minutes, or as infinite time.

In order to set an infinite time it is necessary to go beyond 120 minutes on the display. When these symbols (- - -) appear on the display wait 5 seconds for automatic storage.

ATTENTION: AT INFINITE TIME THE EQUIPMENT CAN BE SWITCHED OFF ONLY MANUALLY.

- **(I) (L) keys**  

The function of these keys is to modify the temperature setting and the incubation time, increasing or decreasing the value. By pressing the **+** key once, with the incubator waiting (B pilot lamp turned on) the set temperature value or the set incubation time is visualized on the display.

From this moment, by pressing the **+** or **--** keys it is possible to modify the value.

After reaching the desired value, wait 5 seconds for automatic acquisition.

If one of these keys is pressed during the cycle, the set values of time or temperature can be read on the display but it is not possible to modify them.

- **(M) Heating pilot lamp**

This lamp is on when resistances are on. After the starting of a cycle, resistances are warm until the set temperature is reached. When the temperature is reached the **HEATING** pilot lamp is turned off and it is alternately turned on, as the equipment is keeping the incubation temperature .

User's instructions

1. Place the material to be treated inside the incubator and close the door (do not open the door during the cycle).
2. Open the two side chimney losses in order to make easier air circulation inside the incubator.
3. Connect the plug into the 220V outlet. Turn the master switch on (A).
4. The "alarm" red pilot lamp is turned on to show the lack of voltage (given by the previous switching off of the machine).
5. Press the "START" push button once, in order to cancel the alarm. The yellow pilot lamp (B) is turned on at this moment, in order to show that the machine is ready for the cycle.
6. Set desired temperature and time of incubation,

Set the temperature : with the temperature pilot lamp on, press the **+** or **--** keys, or Keeping the push button pressed, the increase or the decrease of the set value becomes quick. After reaching the desired value, wait 5 seconds for the acquisition of the set value.

Set the incubation time:

press the switch push button  turning on the timer pilot lamp $\frac{1}{2}$. Repeat the previous operations, in order to visualize and modify the time. Wait 5 seconds for the acquisition of the set value.

7. After the parameters setting, start the incubation cycle by pressing the "START" push button. The yellow pilot lamp (B) is turned off and, if the temperature inside the chamber is lower than the set temperature, the "heating" pilot lamp (M) showing the resistances heating is turned on.
8. The incubator arrives at the set temperature and keeps it for the set time. The temperature measured inside the incubation chamber is visualized on the display (with a $\pm 2^{\circ}\text{C}$ tolerance) or, by switching with the  push button, the remaining time to the end of the cycle.
9. At the end of the set time, the cycle automatically stops and the display shows the writing END, informing that the incubator has finished the cycle.

Important: the cycle is considered valid only if, at the end, the writing END is appearing.

10. At the end of the cycle close the chimney losses in case some material is left inside.
11. If the equipment is not used, turn the master switch (A) off.

ALARM TABLE

| Display/Defect | Cause | Remedy |
|--|-------------------------|---|
| "T.Lo" / inside temperature too low. The minimum heating parameters have not been reached. | 1) lack of voltage | 1) repeat the cycle |
| | 2) excessive load | 2) reduce the load and repeat the cycle |
| | 3) resistances breaking | 3) replace the faulty resistance |
| "T.Hi" / High temperature beyond 110°C | Thermostat breaking | 1) open immediately the door 2) replace the thermostat |
| The equipment is completely stopped | The fuses is broken | Replace the fuses |
| "CnF" | Faulty card | Replace the electronic card |

Cleaning

To carry out a good cleaning use only neutral detergents and rinse with plenty of water.

The gaskets and all plastic parts can be cleaned with a cotton flock soaked in alcohol.

Do not use sponges or abrasive materials which can damage the surface.

To clean do not use detergents, disinfectants or other liquids containing chlorides as they are corrosive for stainless steel; in case of contact please rinse immediately.

Warranty

C.B.M. guarantees its incubator for one (1) years, in case of manufacturing defects or failures of components, not due, according to C.B.M. evaluation, to the operators' wrongful use of the apparatus.

The manufacturing company is released by any responsibilities for damages to people or things due to failure to earth the plant.

The warranty does not cover parts subject to wear and repairs or replacements carried out by personnel not authorised by C.B.M.

IMPORTANT: The manufacturing company is released by any responsibilities for damages to people or things due to failure to earth the plant.

Details, descriptions and pictures are merely indicative.

C.B.M. S.r.l. reserves the right to carry out ameliorative modifications to its own products without any previous notice

Disposal of electric and electronic refuse

Waste from Electrical and Electronic Equipment - Directive 2002/96/EC (WEEE)

The symbol (shown above) that appears on this product means that it should not be disposed of together with normal household refuse but should be treated separately.

Electric and electronic equipment may contain substances that are harmful to the environment and to the health of humans and should therefore be disposed of at special waste disposal sites.

There are harsh legal penalties in force against persons found disposing of such equipment in an irresponsible manner.

When the equipment is no longer of any use, it may taken back to the outlet from which you intend to purchase a replacement.

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