
famos

COMPETENCE IN HIGH TECH SEALING

**Handleiding  
Bedienungshandbuch  
Manual**

**F110**

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## Contenance

0)	Certificate of Conformity and Origin	22
1)	In general	23
	Working area	23
	Installation	24
2)	Operating panel	25
	Menu structure	26
	Explanation of menu structure	26
3)	Preventative maintenance	27
	Technical details	27

© Famos bv, Venlo. The right to make changes is reserved.

# **Certificate of Conformity**

**and**

# **Certificate of Origin**

## **EC declaration of compliance**

In accordance with Appendix 11 sub A of Directive 98/37/EG which dates from 1998, we hereby declare that the listed products conforms with the published specification and complies with the requirements of the following European Directives and Standards and to their current amendments:

- Machines guideline 98/37/EG
- European standard EN 60335-1:2002 + A:112004 + A1:2004  
EN 60335-2-45:2003  
EN 50366:2003
- EMC guideline (Richtlinien 89/336/EEG)
- Low voltage guideline 2006/95/EG

## **Product**

Type

F110 serie

from 01-01-2007

## **Certificate of origin**

Herewith it is certified that the listed equipment has been manufactured in a member state (The Netherlands) of the European Community.

Famos bv

T.J.E. van Gardingen  
General Director

## **1) In general**

All Famos rotary heat sealers conform to the Dutch and German Directives (R3280 en DIN 58.953 Teil 7) that are drawn up for sealing packing material as indicated below. All Famos F110 heat sealers can therefore be used for packing medical equipment.

The Famos heat sealers are **NOT** medical products.

It is not allowed to use the sealers in the direct presence of patients (for example, in operation theatres). The minimum distance between the heat sealing machines and patients should be 1,5m.

## **Working area**

The Famos F110 rotary heat sealer is suitable for sealing following materials:

- Laminate packing material
- Sterilisation paper packingmaterial fitted with top seal
- Tyvek (if sealer is designed as a Tyvek sealer)

For the correct temperature setting of the different packing materials we refer to the manufacturer of the materials.

The Famos F110 rotary heat sealer is **not** suitable for sealing the following material:

- Polyethylen foil
- Soft PVC foil
- Hard PVC foil
- Polyamide foil
- Polypropylen foil

(For above materials Famos bv can offer an impulse sealer).

## **Installation**

In order to guarantee a correct functioning, we advise to install the sealer at a spot:

- where the machine has enough space for ventilation.
- where it is not exposed to direct sunlight or some other heatsource.
- where it is not subject to vibrations or shocks.
- that is dustfree.

Check if the on/off switch is off.

### REMARK


Connect the mains lead to the sealer first and than to an earthed socket with the correct Mains voltage (according to type label at the rearside of the sealer).

Switch on the on/off switch.

During the heating-up the machine will mention in the display that the temperature is too low. From the off-situation onwards it takes  $\pm 5$  minutes before the pre-set temperature (at 180°) is reached.

Two beeps indicate that the pre-set temperature is reached. The display will show the stand-by screen.


Feeding packages is possible now.

## 2) Operating

### Operating panel

### Meaning of the keys

- |   |                       |   |
|---|-----------------------|---|
| A | Menu keys             | Selecting a menu or option depicted above the key.  |
| B | Cursor keys           | Moving the cursor in the direction of the arrow.  |
| C | Programming key       | Access to the main menu.  |
| D | Transport key         | Transporting packages to the left.  |
| E | Alphanumeric keyboard | 0... 9, A...Z, ( [space] press 4x key 9), ( [.] press 2x key 0), ( [-] press 3x key 0), ( [:] press 4x key 0) |
| F | Enter key             | Store data / Display previous menu.   |
| G | Counter reset key     | Resetting the counter to zero.  |
| H | Delete key            | Deleting a character on the cursorposition (alphanumeric entrance)  |
| I | Insert key            | Inserting a characterspace on the cursorposition (alphanumeric entrance)                                      |



### Entering only numerical information using the keyboard

By pressing a number key the relevant number is entered and the cursor jumps one position to the right. The number can be stored by pressing the enter key

### Selecting one option

By pressing the menu key once under the current value, this value changes into the other value (Yes $\diamond$ No, Wide $\diamond$ Narrow ) or the relevant setting can be selected (DMY, YMD, MDY).

**Menu structure**

**Explanation :**

**A**  
Sys                    Setting the system data. ( > = scroll menu )

**B**  
Temp                  Temperature setting.  
Dat/Time              Setting date and time information.  
Misc.                   Various settings.

**C**  
Time                   Setting the time (24 hour system).  
Date                    Setting the date.  
Structure               The way in which the date is indicated.  
                              DMY=day,month,year |YMD=year/month/day |MDY=month,day,year  
Counter                Display of the total number of packages sealed. The sealmachine gives 4 beeps  
                              after the number of packages – entered here – has been sealed.  
Language                Selecting the menu language ( > = scroll menu )

**Preventative maintenance**

In order to guarantee an optimal functioning of the rotating sealer and in order to prevent corrective maintenance, we advise to have the machine checked and cleaned at least twice a year. During this inspection the following parts are checked on wear and if necessary replaced.

- 12.210.367 Transportbelts
- 30.200.039 Teflon strip on heating elements
- 30.200.041 Teflon strip on guiding bars

For the article numbers and descriptions of all other parts we refer to your dealer and/or to the technical documentation.

### **Technical details**

- Temperature range : 100 - 199 °C
- Feeder speed : 8 m/min
- Seal width : 15 mm.
- Self adjusting pressure : 10 Kg (Intern mechanically adjustable: 1 – 12 Kg)
- Fuses : See type label at the rearside of the sealer
- Capacity : max. 365 W, average 165 W
- Main voltage : See type label at the rearside of the sealer
- Dimensions : 550x320x210 mm
- Weight : 18,5 Kg

#### Warning

To avoid the risk of electric shock do not open the casing. Repairs should be done by a specialist.