

Hangzhou Zhengda Medical Co., Ltd

Producing the article protects to fix single

Ltd. model:

Number of the machine:

Buy date:

User's name: _____

User's address _____

Telephone number _____

Invoice number:

Maintain the record						
Date	The reason of Maintain	The contain of Maintain	The cost of Maintain	Maintain signs	User's signs	Remark

Precautions

1. The products free guarantee period that you bought is half a year from day when buy.
2. If the products that you bought have any doubtful point while using, guidance and service.
3. If because users or exceeds guarantee period, our company. Collects the fee of material only according to the circumstances.
4. Must return this card while main taining, so please keep this singly.

Address:No.279 Shiqiao Road, Hangzhou, China
 Postcode:310022
 Maintenance Phone:0571-85190105

Operation method and precautions

1, Insert the plug in the outlet, hold the gypsum tightly to saw all alone. Notice that the frequency of gypsum saw is more than 133hz, so before using, don't close the switch of gypsum saw to avoid dangerous.

2, Hold the gypsum with a firm grip and open the switch, if there is no abnormality, put saw blade on the surface which need cutting, then the vertical gypsum is cut down, if the gypsum material is the hard polyester material, you can move the saw parallelly to the gypsum surface. If not sure, don't cut entirely (although this gypsum has the protect function for soft tissue).

3, vertically cut down the gypsum, lift the gypsum and move forward, then cut the next section of gypsum saw surface.

4, Repeat the above-mentioned course, until finish the task.

5, If the cut section of saw is worn and torn, can unscrew the set bolt with the special-purpose spanner, rotate an angle to tighten it again.

6, It is the same to the method of changing the saw blade as above-mentioned..

7, If the frequency of output drops while using, such as the noise is great, should generally stop continuing

8, Check the cable regularly while using to avoid using torn cable.

Disinfect method and precautions

This Electric Plaster Cutting Saw is strictly forbidden using the high temperature, high-pressure method to disinfect, can use soft cotton cloth and neutral soap liquid to disinfect it. Forbid using the organic solvent and rigid material. Don't let the liquid into the gypsum saw while cleaning.

Simple failure elimination and maintenance

If this Electric Plaster Cutting Saw can't run, please check first whether the carbon brush is worn and torn, in this way, the ones that should adopt our company's special using carbon brush for changing to prevent the armature from wearing and tearing abnormally. If there are relevant problems, please contact us.

If there is abnormal failure in operating, it may lack lubricating grease which can bear the high-temperature, can get in touch with us.

Should notice to dry it with hair-dryer after using, in order to protect the electrical machinery, lengthen its span life. If there are other problems, please getting in touch with us.

Commitment for after service

If the product has any quality problem in one year, we provide the service of change and repair, and we also provide long-term repairing service. And only charge the material expenses in guarantee period.

Hangzhou Zhengda Medical Co.,Ltd

Address: No.279 Shiqiao Road, Hangzhou, China

Tel: 0571-85190105, 85096839

Postal code: 310022

E-mail: nwsales@zd-

medicail.com

Website: www.zd-medical.com


Technique feature of product

YTJ Electric Plaster Cutting saw is developed for the second generation products of our company; it is the newly-designed medical equipment products totally different from the first generation. Compared to the first generation of products, this new product considers the life-span of new material (such as the polyester) and medical use instruments. It is revolutionary (been improved from 20HZ to more than 133 HZ in the technical design of operation frequency. it is more swift to cut. Distance slice of the new instrument has adopted the new material; each of its technical indicators is higher than the first generation.

The main capability parameters are as follows:

The maximum unload Output rev :	>10000R.P.M
Swing angle:	$\geq 4^\circ$
Rated input voltage :	220Vac \pm 22Vac
Rated power:	<125W
Noise:	<85dB (When being unloaded)
Weight:	<1.5Kg
Length:	<250mm

Graphic presentation of electric gypsum saw:



Y&Z

Please read the manual before using

OWNER'S MANUAL

Electric Plaster Cutting Saw

Type

YTJ

Hangzhou Zhengda Medical Co., Ltd

This publication is printed



On 100% recycled paper