

CBC Kimetec Tourniquets

Latex Free

Designed with patient comfort in mind, this tourniquet offers a velvety-soft touch. It's easy-to-use and features the popular one-hand operation, which prevents pinching, catching, and jamming. Plus, you'll love the quick release buckle that allows for instant removal. It's a great way to put the patient at ease and make your job easier!

Features

- Skin-friendly, flexible and soft ribbon 25mm wide
- Sturdy safety snap closure
- One-hand operation
- Available with quick release lock
- Metal Free and Latex Free
- No squeezing or pinching of the skin
- CE, Made in Germany



Code	Description
K5020	TOURNIQUET MULTI-COLOURED DESIGNER
K5085	TOURNIQUET DINO DESIGN
K5087	TOURNIQUET KIDS DESIGN
K5086	TOURNIQUET CLOWN DESIGN
K5660	TOURNIQUET INFANTS BLUE (SIZE: L32CM X W1.0CM)
K5670	TOURNIQUET INFANTS RED (SIZE: L32CM X W1.0CM)



Code	Description
K5030	TOURNIQUET PINK
K5040	TOURNIQUET MINT
K5050	TOURNIQUET GREY
K5060	TOURNIQUET BLUE
K5070	TOURNIQUET RED
K5080	TOURNIQUET ORANGE
K5090	TOURNIQUET GREEN



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: K50_Tourniquets-CBC-Kimetec_v4-2

First Prepared: 13/11/14 Eevien Tan
 Release Approval: 17/11/14 Lee Carter
 Version Control: V4.2 Page 1 | 3

CBC Kimetec Tourniquets

Latex Free

Material					
Casing	Polycarbonate (PC)				
All Other Casing Parts	Polyoxymethylene (POM)				
Ribbon	Polyester and Lycra®, phthalate free				
Colours of Ribbon	The colours do NOT contain the following chemical substances which are prohibited in textiles: 4-Aminodiphenyl, Benzidine, 4-Chlor-o-Toluidine, 2-Naphthylamine, o-Aminoazotoluene, p-Chloraniline, 2-Amino-4-nitrotoluene, 2,4-Diaminoanisole, 4,4-Diaminodiphenylmethane, p-Cresidine, 4,4-Methylen-bis-(2-chloraniline), 4-4-Oxydianiline, 4,4-Thiodianiline, o-Toluidine, 2,4-Toluenediamine, 2,4,5-Trimethylaniline				
Specifications to Maintain Hygiene					
Rating	The tourniquet can be rated as uncritical medical product, provided that it is applied on intact skin. Do not apply on pathologically altered skin.				
Cleaning Procedure					
Washing	With minor staining, the tourniquet can be washed by hand with domestic laundry detergents at up to 60°C. It is recommended to use hand wash detergents which offer a good cleansing effect at a temperature of 40°C. Please observe the respective manufacturer's instructions of use. Please rinse the tourniquet at least twice with fresh, lukewarm water. Avoid scented detergents. Due to the variety of detergents available in the market, we cannot recommend a specific product.				
Note	If the tourniquet is contaminated with blood, secretion, or other body fluids, please dispose of it immediately in accordance with all pertinent regulations on the disposal of biohazardous waste. In this case, it is not permissible to clean and reuse the tourniquet.				
Procedure of Disinfection					
Physical Procedure	Boiling in water/steam sterilisatoin				
Note	On the part of the manufacturer, the product is not suited for these disinfection procedures.				
Chemical Procedure	For our tourniquet, we recommend a disinfection with a liquid disinfectant. Use CIDEX OPA by Medos International Sarl as disinfectant. Place the tourniquet in a container. Do not use a metallic container. Do not coil up the tourniquet. Pour as much disinfectant solution as required into the container to completely cover the tourniquet and leave it for at least 5 minutes in the solution. The room temperature shall be at least 20°C. After the contact time, rinse the tourniquet with plenty of fresh water for at least 3 rinsing rounds each with fresh water. Leave the tourniquet in the rinse water for at least 1 minute. You can hang the tourniquet up for drying. Make sur that dripping water does not do any damage. For more information on safe use, disposal, and storage life of the disinfectant, refer to the corresponding instructions of use.				
Manufactured and distributed by Cidex Opa	<table border="0"> <tr> <td>Manufacturing:</td> <td>Distribution in Germany:</td> </tr> <tr> <td>Medos International Sarl Chemin-Blanc 38 2400 Le Locole, Switzerland</td> <td>Johnson+Johnson Medical GmbH Oststraße 1 22844 Norderstedt</td> </tr> </table>	Manufacturing:	Distribution in Germany:	Medos International Sarl Chemin-Blanc 38 2400 Le Locole, Switzerland	Johnson+Johnson Medical GmbH Oststraße 1 22844 Norderstedt
Manufacturing:	Distribution in Germany:				
Medos International Sarl Chemin-Blanc 38 2400 Le Locole, Switzerland	Johnson+Johnson Medical GmbH Oststraße 1 22844 Norderstedt				



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: K50_Tourniquets-CBC-Kimetec_v4-2

First Prepared: 13/11/14 Eevien Tan
Release Approval: 17/11/14 Lee Carter
Version Control: V4.2 Page 2 | 3


Disclaimer: All images are included for reference only. Features and specifications are subject to change without prior notice.
Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes. E & OE
Visit our website at: <http://www.lrinstruments.com.au>

CBC Kimetec Tourniquets


Latex Free

Instructions For Use


To ensure a safe application of the product, we recommend that it should be used by medical staff only.




Step 1:
Close




Step 2:
Tie Up



Step 3:
Slow Release



Step 4:
Hold Down, Then Open



Step 5:
Roll Up



Kimetec GmbH
Gerlinger Straße 36-38
71254 Ditzingen Germany
www.kimetec.de



Conformity mark according to EU Council Directive 93/42/EEC
GB2010001M
Stand/As of 2014-04

Cleaning: The product can be washed with gentle fabric detergents.

Disinfection: Detailed information is available on request.

NOTE: The tourniquet can be autoclaved 12 times at 121°C and 1.2 bar for 20 minutes. 134°C is not possible.



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: K50_Tourniquets-CBC-Kimetec_v4-2

First Prepared: 13/11/14 Eevien Tan
Release Approval: 17/11/14 Lee Carter
Version Control: V4.2 Page 3 | 3

Disclaimer: All images are included for reference only. Features and specifications are subject to change without prior notice.
Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes. E & OE
Visit our website at: <http://www.lrinstruments.com.au>