

FAMOS Competence in CSSD

Hardware & Sterilisation Consumables

Product Information

A single company with a reputation for Technical Excellence and Product Reliability, committed to CSSD'S Manufacturing both Hardware and Sterilisation Consumables

CSSD Hardware Manufactured in Holland

- Instrument Baskets
- Workbenches
- Heat Sealers



Lot Control
Standard Cycles
B & S Cycles
EN867 Compliant
Class 1 Medical Device TGA Incl



Steam integrator - Prion
EN867 Class D
ISO1140-1 Class 6
Class 1 Medical Device TGA Incl

Sterilisation Consumables Manufactured in Germany

- Integrators
- Bowie Dicks
- Lot Controls



Autoclave Tape
19mm - 25mm Widths
Extra Adhesive Available
EN867 - Class A Compliant
Class 1 Medical Device TGA Incl



Class 6 Steam Integrator
EN867 Class D
ISO1140-1 Class 6
Class 1 Medical Device TGA Incl

What is your Sterilisation Assurance Plan?
No two Sterilisation Quality Assurance Systems are the same!
Talk to us and we will work with you to provide for your particular needs.

	Manufactures both CSSD Hardware & Consumables both CSSD Hardware & Consumables	Bowie dick Air Removal Tests	Specialist Lot Controls for Hospital Sterilisers	Specialist Lot Controls for Class B & S Cycles. Ideal for Dental & GP Clinics	Adhesive Labels for Plasma Sterilisation	Produces Customer Specific Labels & SAP (Sterility Assurance Products)	Lot Control for PRION Cycles (134 degrees, 18 minutes)	Integrators for PRION Cycles (134 degrees, 18 minutes)
FAMOS	YES	YES	YES	YES	YES	YES	YES	YES

Wherever Sterilisation Assurance is needed, including:

- | | |
|------------------|---------------------|
| Hospitals | GP Surgeries |
| Laboratories | Clinics |
| Dental Surgeries | eterinary Surgeries |

