



BOWIE-DICK TEST PACK STERILISATION MONITORING

Bowie Dick Test Pack for pre-vacuum hospital steam sterilisers (porous loads only)

The benefits of a Traditional Single Use Bowie Dick Test Pack, with the addition of easy and unique Documentation Boxes incorporated on to the sheet and a 5 year shelf life from date of manufacture. The product conforms to EN 867-4 and ISO 11140-4 and is available in packs of 20.

The Bowie-Dick Test has been a mandatory test for steam sterilisers for quite some time and is a functionality test that is required to be completed before one can release the steam steriliser for use for that day. Based upon the load configuration the selection in the type of Bowie-Dick test can be made. Famos manufacturers both the Bowie-Dick test pack (for porous loads) and Bowie-Dick test set based upon a helix device (for hollow loads).

- Indicator complies to ISO 11140 Class 2
- Good and easy documentation by attached data input areas
- Easy to interpret colour change from yellow to black
- Laminated indicator prevents any interaction with ink
- Paper is high quality re-cycled paper
- Ink is free of toxic components
- Non reversible indicator



CODE	DESCRIPTION
F15.100.069	BOWIE-DICK TEST PACK FOR POROUS LOADS (PK 20)

E&OE First Published: 14 April 2009
Last Updated: 22 February 2011
File Name: info_F15.100.069_BOWIE-DICK-TEST-PACK

