



BOWIE-DICK HELIX TEST STERILISATION MONITORING

Bowie-Dick Helix Test for:

- (1) Pre-Vacuum Hospital Steam Sterilisers (hollow load only)
- (2) Bench-top Steam Sterilisers Class B

Small bench top steam sterilisers cannot be tested with a Bowie-Dick Test pack. Therefore EN 867-Part 5 and ISO 11140 defines a Helix PCD based Steam Penetration Test as an alternative for the Bowie-Dick Test.

Famos designed a red lid for the Helix PCD device and marked the chemical indicators with red, not to let any mix-up with lot control PCD and indicators.

- Helix PCD complies to EN867 part 5 and ISO 11140
- Indicators comply to ISO 11140 Class 2
- Easy to use and document
- Non reversible indicator
- Accurate and reproducible results
- Best interpretation of color change from yellow to black
- Suitable to be used in conjunction with the SASCI
- Ink is free of any toxic component
- Cost effective test
- Environmentally friendly



CODE	DESCRIPTION
F15.100.097	BOWIE-DICK HELIX TEST FOR HOLLOW LOADS (PK 100)
F15.100.098	BOWIE-DICK HELIX TEST FOR HOLLOW LOADS (PK 250)

E&OE First Published: 14 April 2009
 Last Updated: 22 February 2011
 File Name: info_F15.100.097_BOWIE-DICK-HELIX-TEST

