

do YOU
conform
to 4187?

AS/NZS 4187:2003
Australian/New Zealand Standard™
Cleaning, disinfecting and sterilizing
reusable medical and surgical
instruments and equipment, and
maintenance of associated
environments in health care facilities



do YOU use a
quiver
with a
removeable
base?



do YOU use the
correct
instrument
pin?



Quivers

Protect everyone from inadvertent fires, use a Quiver!

- Removable Base
- Autoclavable
- Latex Free

Detachable for easy cleaning and sterilization. Available in 16, 20, 25, 30, 35 and 40cm lengths to safely manage open and laparoscopic procedures.



Instrument Pins

AS/NZS 4187:2003 states:

“Instruments must be MAINTAINED open during all stages of reprocessing”.

Using an instrument pin, effectively reduces handling time, reprocessing time and ensures you meet 4187.



Backhaus Towel Clips

Available in varying sizes, ranging from 5cm in length to 19.5cm in length.

- Stainless Steel
- Secure Ratchet Locking
- Backhaus Style



Product Information




BASKETS - PERFORATED SHEET

Code	Description
CBM75	PERFORATED SHEET BASKET STAINLESS STEEL 540 X 240 X 45
CBM75/1	PERFORATED SHEET BASKET STAINLESS STEEL 540 X 240 X 50
CBM75/2	PERFORATED SHEET BASKET STAINLESS STEEL 540 X 240 X 70
CBM75/3	PERFORATED SHEET BASKET STAINLESS STEEL 540 X 240 X 90
CBM75/4	PERFORATED SHEET BASKET STAINLESS STEEL 540 X 240 X 41
CBM75/5	PERFORATED SHEET BASKET STAINLESS STEEL 540 X 240 X 180
CBM78	PERFORATED SHEET BASKET STAINLESS STEEL 440 X 340 X 150
CBM78/1	PERFORATED SHEET BASKET STAINLESS STEEL 440 X 340 X 70
CBM78/2	PERFORATED SHEET BASKET STAINLESS STEEL 440 X 340 X 200
CBM79/2	PERFORATED SHEET BASKET STAINLESS STEEL 650 X 240 X 70



BASKETS - WIRE MESH

Code	Description
CBM84	WIRE MESH BASKET STAINLESS STEEL 480 X 240 X 50
CBM84/1	WIRE MESH BASKET STAINLESS STEEL 480 X 240 X 70
CBM84/2	WIRE MESH BASKET STAINLESS STEEL 480 X 240 X 50
CBM85	WIRE MESH BASKET STAINLESS STEEL 240 X 240 X 50
CBM85/1	WIRE MESH BASKET STAINLESS STEEL 240 X 240 X 70
CBM85/2	WIRE MESH BASKET STAINLESS STEEL 240 X 240 X 90
CBM86	WIRE MESH BASKET STAINLESS STEEL 540 X 240 X 50
CBM87	WIRE MESH BASKET STAINLESS STEEL 540 X 240 X 90
CBM88	WIRE MESH BASKET STAINLESS STEEL 210 X 145 X 50
CBM89	WIRE MESH BASKET STAINLESS STEEL 300 X 240 X 50
CBM90	WIRE MESH BASKET STAINLESS STEEL 540 X 240 X 70



BASKETS - MINI

- Designed to accommodate those small items that get lost in a tray and for transportation throughout washing process
- Lids overlap to accommodate easy opening with a gloved hand

Code	Description
CBM63T	CBM MINI BASKET 170X50X33MM INCL SILICONE MAT
CBM67S	CBM MINI BASKET 170X135X35MM INCL SILICONE MAT
CBM77T	CBM MINI BASKET 170X270X35MM INCL SILICONE MAT