



ACCESSORIES FOR MEDICAL & SURGICAL INSTRUMENTS

do YOU conform to 4187?

AS/NZS 4187:2003
Australian/New Zealand Standard™

Cleaning, disinfecting and sterilizing reusable medical and surgical instruments and equipment, and maintenance of associated environments in health care facilities

Standards Australia STANDARDS NEW ZEALAND

do YOU use the correct instrument pin?

do YOU use a quiver with a removeable base?

INSTRUMENT PINS

Instrument pins or sterilising pins are designed to hold surgical instruments during the sterilisation process. Surgical instruments are placed in the open position on the pins keeping the instruments in the most advantageous position for proper sterilisation.

AS/NZS 4187:2003 states: "Instruments must be MAINTAINED open during all stages of reprocessing". Using an instrument pin, effectively reduces handling time, reprocessing time and ensures you meet 4187.



Instruments are placed in the open position on Sterilising Pin (Instruments not included)

CODE	DESCRIPTION
HRA-002-004-PMS	CHIFA STERILISING PIN 102 X 63MM
HRA-002-006-PMS	CHIFA STERILISING PIN 152 X 63MM
HRA-002-008-PMS	CHIFA STERILISING PIN 203 X 63MM
HRA-002-010-PMS	CHIFA STERILISING PIN 254 X 63MM
HRA-002-012-PMS	CHIFA STERILISING PIN 300 X 63MM
LIH100X65	LIBERTY STERILISING PIN 100 X 65MM
LIH150X65	LIBERTY STERILISING PIN 150 X 65MM
LIH200X65	LIBERTY STERILISING PIN 200 X 65MM

E&OE First Published: 1 August 2010
Last Modified: 3 November 2010
File Name: info_LIH_INSTRUMENT-PINS

