

# MR-100 Manual Resuscitators



## Adult Resuscitator

- Autoclavable to 170°C
- Specifications: MR-100 Adult
- See-through carry case with instruction manual

Material	Polysulfone Silicone Rubber Polyvinyl chloride
Connectors	15/22mm ISO patient connector 24mm ID silicone neck 26mm ID (to silicone bag inlet) 25mm OD (to oxygen reservoir) 6mm oxygen inlet
Inspiratory/ Expiratory Resistance	2.5cm H2O/3.1cm H2O at 50 lpm flow
Operating Conditions	-18°C to 50°C
Reservoir Volume	2500ml
Maximum Oxygen Concentration	99% (with oxygen supply)
Pressure Relief	60 ± 10cm H2O lock type override mechanism
Maximum/Stroke Volume	1500/1350ml
Maximum Breath Cycle Rate	45 bpm (600ml)



Code	Description	Mask	Reservoir	Pop-Off
1068	MR-100 ADULT RESUSCITATOR	5145	2500ML	60CM H2O
0043	ADULT SEE-THROUGH CARRY CASE			



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107  
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134  
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au  
File Name: MR-100\_Manual-Resuscitators\_v2-0

First Prepared: 05/02/15 Eevien Tan  
Release Approval: 05/02/15 Lee Carter  
Version Control: V2.0 Page 1 | 2

# MR-100 Manual Resuscitators



## Child Resuscitator

- Autoclavable to 170°C
- Specifications: MR-100 Child
- See-through carry case with instruction manual

Material	Polysulfone Silicone Rubber Polyvinyl chloride
Connectors	15/22mm ISO patient connector 24mm ID silicone neck 26mm ID (to silicone bag inlet) 25mm OD (to oxygen reservoir) 6mm oxygen inlet
Inspiratory/ Expiratory Resistance	2.5cm H2O/3.1cm H2O at 50 lpm flow
Operating Conditions	-18°C to 50°C
Reservoir Volume	2500ml
Maximum Oxygen Concentration	99% (with oxygen supply)
Pressure Relief	40 ± 5cm H2O
Maximum/Stroke Volume	550/350ml
Maximum Breath Cycle Rate	105 bpm (70ml)



Code	Description	Mask	Reservoir	Pop-Off
1086	MR-100 CHILD RESUSCITATOR	-	-	40CM H2O
0041	CHILD SEE-THROUGH CARRY CASE			



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107  
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134  
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au  
File Name: MR-100\_Manual-Resuscitators\_v2-0

First Prepared: 05/02/15 Eevien Tan  
Release Approval: 05/02/15 Lee Carter  
Version Control: V2.0 Page 2 | 2