

CPR Super Pocket Resuscitator



Award-winning and patented product. Combines extension tube, one-way valve, foldable face mask and oxygen inlet in a small and compact clam-shell carrying case. Best tool for mouth to mask resuscitation. Minimises cross-contamination during training or CPR practice. High efficiency one-way valve with swivel feature and standard 15M connection. Easy to access - ergonomic mask with soft air-cushion for comfort and better seal and flexible extension tube for better adjustment. Easy to carry in an all-in-one compact storage case.

Specifications

Connections	15mm OD
Inspiratory Resistance	4.5cm H2O
Expiratory Resistance	3.5cm H2O
Material	PVC - Single patient use
Dead Space	125ml
Oxygen Concentration	40-80% (depending on flow and ventilation cycles)
Measurements (case)	132 x 107 x 46 mm
Weight	100g



Features

- Bacterial air filter with leakage free air cushion mask
- Silicone duck bill one way valve, Oxygen port
- Disposable: Single Use Only
- Detachable mouthpiece - mask can then be connected to other resuscitation devices
- Adult size only, also suitable for large child
- Mask shell folds for compact storage within the compact carry case
- Custom labelling - bulk quantities apply
- Sold individually
- Packaged in boxes of 4



Code	Description
1325	CPR SUPER POCKET RESUSCITATOR WITH FILTER

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: 1325_Pocket-Resuscitator_v2-0



First Prepared: 22/01/15 Eevien Tan
 Release Approval: 22/01/15 Lee Carter
 Version Control: V2.0 Page 1 | 1