

# Manual Silicone Resuscitators



## Adult Resuscitator

Specifications: RESCU-2 Adult

- Autoclavable to 121°C
- Includes instruction manual

Material	Polycarbonate, Silicone Rubber Polyvinyl Chloride
Inspiratory / Expiratory Resistance	<5cm H2O at 50 lpm flow
Reservoir Volume	2500ml
Max. Oxygen Concentration	99% (with oxygen supply)
Operating Conditions	-18°C to 50°C
Connectors	15/22 mm ISO patient connector 24 mmID silicone neck 26 mmID (to silicone bag inlet) 32 mmID (to combined oxygen valve & reservoir)
Pressure Relief	60 ± 10cm H2O Lock type override mechanism
Maximum/Stroke Volume	1500/1350ml
Max Breath Cycle Rate	45 bpm
Mask	Up to 25 autoclave cycles
Reservoir	Single use, disposable
Expiry	3 Years from date of manufacture
Optional Extra	See-through carry case



Code	Description	Mask	Reservoir	Pop-Off
1405	SILICONE ADULT RESUSCITATOR	5135	2500ML	60CM H2O
0043	ADULT SEE-THROUGH CARRY CASE			
0041	CHILD SEE-THROUGH CARRY CASE			



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107  
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134  
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au  
File Name: 140\_Manual-Resuscitators\_v4-2

First Prepared: 22/01/15 Eevien Tan  
Release Approval: 22/01/15 Lee Carter  
Version Control: V4.2 Page 1 | 2

# Manual Silicone Resuscitators



## Child Resuscitator

Specifications: RESCU-2 Child

- Autoclavable to 134°C
- Includes instruction manual

Material	Polycarbonate, Silicone Rubber Polyvinyl Chloride
Inspiratory / Expiratory Resistance	<5cm H2O at 50 lpm flow
Reservoir Volume	2500ml
Max. Oxygen Concentration	99% (with oxygen supply)
Operating Conditions	-18°C to 50°C
Connectors	15/22 mm ISO patient connector 24 mmID silicone neck 26 mmID (to silicone bag inlet) 32 mmID (to combined oxygen valve & reservoir)
Pressure Relief	40 ± 5cm H2O
Maximum/Stroke volume	1500/1350ml
Breath Cycle Rate	105 bpm
Mask	Up to 25 autoclave cycles
Reservoir	Single use, disposable
Expiry	3 Years from date of manufacture
Optional Extra	See-through carry case



Code	Description	Mask	Reservoir	Pop-Off
1406	SILICONE CHILD RESUSCITATOR	5122	2500ML	40CM H2O



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107  
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134  
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au  
File Name: 140\_Manual-Resuscitators\_v4-2

First Prepared: 22/01/15 Eevien Tan  
Release Approval: 22/01/15 Lee Carter  
Version Control: V4.2 Page 2 | 2