

Oxygen Masks

- Soft, transparent, odour free vinyl
- Aluminium fixer and elastic strip to ensure the best fit and user comfort



Code	Description
3644	OXYGEN MASK CHILD & 2M O2 TUBING
3645	OXYGEN MASK ADULT & 2M O2 TUBING

Hi-Oxygen Concentration Masks

- Soft, transparent, odour free vinyl
- Aluminium fixer and elastic strip to ensure the best fit and user comfort
- One way valve and reservoir to increase the efficiency of oxygen concentration



Code	Description
3674	HIGH OXYGEN MASK CHILD WITH 2M TUBING AND RESERVOIR
3675	HIGH OXYGEN MASK ADULT WITH 2M TUBING AND RESERVOIR

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: 36_Oxygen-Masks_v3-1



First Prepared: 22/01/15 Eevien Tan
Release Approval: 22/01/15 Lee Carter
Version Control: V3.1 Page 1 | 1