

Bacterial Filter With O2 Port

Features

- Designed for a better inlet air quality
- Designed to prevent cross contamination and protection of the operator
- Efficiency: 99% for 0.3 micron (3M's material)
- Connector 22 mm ID and 22/15 mm



Code	Description
39016	DISPOSABLE BACTERIAL FILTER WITH LEUR PORT
39014	HEPA RADIAL PLEATED FILTER
70015	BACTERIAL FILTER WITH CLEAR HOUSING (NOT PICTURED)

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: 3901X_Bacterial-Filters_v3-1



First Prepared: 22/01/15 Eevien Tan
 Release Approval: 22/01/15 Lee Carter
 Version Control: V3.1 Page 1 | 1