

# Littmann Classic II SE Stethoscope

Featuring a two-sided chest-piece, this classic stethoscope offers high acoustic sensitivity for superior performance. A tunable diaphragm on one side augments the traditional bell function of the opposite side.

The Littmann Classic II S.E. stethoscope offers a traditional bell combined with a tunable diaphragm that allows the user to conveniently alternate between low- and high-frequency sounds without turning over the chest-piece.

## Features

- Traditional bell combined with a tunable diaphragm that allows the user to conveniently alternate between low and high-frequency sounds without turning over the chest-piece
- Excellent acoustic seal and comfortable fit with patented 3M™ Littmann® Snap Tight Soft-Sealing Ear-tips
- Comfortably angled, anatomically correct headset
- Patient-friendly, non-chill rim and diaphragm
- Single-lumen tubing is available in a variety of popular colours
- Three-year warranty



Ear Tips



Turnable Diaphragm



Code	Description
3M2201	LITTMANN CLASSIC II SE STETHOSCOPE, BLACK
3M37809P	LITTMANN EAR TIP SOFT SNAP, BLACK SMALL PAIR
3M37811P	LITTMANN EAR TIP SOFT SNAP, BLACK LARGE PAIR
3M36556S	LITTMANN DIAPHRAGM/RIM SET TURNABLE, BLACK SINGLE

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107  
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134  
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au  
 File Name: 3M2201\_Littman-Classic\_v3-1



First Prepared: 21/01/15 Eevien Tan  
 Release Approval: 21/01/15 Lee Carter  
 Version Control: V3.1 Page 1 | 1