

Aerosol Therapy



	3736 / 3737	3569 / 3570
Principle of Nebulisation	Venturi	Jet
Maximum Capacity	30ml because the higher capacity is better for long term application or lower dosage	9ml a proper capacity for less than 15 minute treatment
Dead Space (residual volume)	1.5ml	0.5ml big saving for expensive drugs
Particle Size	0.5-5 micron	0.5-5 micron
Deposition (medicine delivery)	0.2ml / minute (7 lpm)	0.2ml / minute (6lpm) higher efficiency)
Flow Range	4-12 lpm	4-12 lpm
Connector (input)	6mm	6mm
Connector (output)	19mm	19mmID & 22mm OD
Position for Nebulisation	Vertical or limited angle	Any position with proper medicine container
Bottle	Transparent milky, leakage free design	Transparent clear, leakage free design



Venturi Nebulisation



Jet Nebulisation



Code	Description
3569	NEBEASY JET NEBULISER MASK & TUBING ADULT
3570	NEBEASY JET NEBULISER MASK & TUBING CHILD
3736	VENTURI NEBULISER MASK & TUBING CHILD
3737	VENTURI NEBULISER MASK & TUBING ADULT
3732	AEROSOL MASK CHILD
3733	AEROSOL MASK ADULT
35170	VENTURI NEBULISER BOWL ONLY
3551	JET NEBULISER BOWL ONLY
31310	OXYGEN TUBING 2 METRES
3002	MOUTH PIECE ONLY

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: Aerosol-Therapy_v2-0



First Prepared: 05/02/15 Eevien Tan
 Release Approval: 05/02/15 Lee Carter
 Version Control: V2.0 Page 1 | 3

Aerosol Therapy



Aerosol Masks | Oxygen

- Soft, transparent, odour free vinyl
- Aluminium fixer and elastic strip to ensure the best fit and user comfort
- Non rebreathing



Code	Description
3642	OXYGEN MASK CHILD WITHOUT TUBING (NOT PICTURED)
3643	OXYGEN MASK ADULT WITHOUT TUBING (NOT PICTURED)
3644	OXYGEN MASK CHILD WITH 2M O2 TUBING
3645	OXYGEN MASK ADULT WITH 2M O2 TUBING



Aerosol Masks | Hi-Oxygen Concentration



- Soft, transparent, odour free vinyl
- Aluminium fixer and elastic strip to ensure the best fit and user comfort
- One way valve and reservoir to increase the efficiency of oxygen concentration
- Non rebreathing



Code	Description
3674	HIGH OXYGEN MASK CHILD WITH 2M TUBING & RESERVOIR
3675	HIGH OXYGEN MASK ADULT WITH 2M TUBING & RESERVOIR

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: Aerosol-Therapy_v2-0



First Prepared: 05/02/15 Eevien Tan
 Release Approval: 05/02/15 Lee Carter
 Version Control: V2.0 Page 2 | 3

Aerosol Therapy



Oxygen Cannula & Tubing

- Soft, transparent, odour free vinyl to ensure the best fit and user comfort
- Flexible oxygen tubing prevents occlusion when tubing is kinked



Code	Description
3231	NASAL CANNULA INFANT WITH 2M TUBING
3232	NASAL CANNULA CHILD WITH 2M TUBING
3233	NASAL CANNULA ADULT WITH 2M TUBING

Aerosol Therapy Accessories



Code	Description
35170	VENTURI NEBULISER BOWL ONLY
3551	JET NEBULISER BOWL ONLY
75047	T-PIECE ONLY
3002	MOUTH PIECE ONLY
31310	OXYGEN TUBING 2 METRES
31610	OXYGEN TUBING 10 METRES



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: Aerosol-Therapy_v2-0

First Prepared: 05/02/15 Eevien Tan
 Release Approval: 05/02/15 Lee Carter
 Version Control: V2.0 Page 3 | 3