

Basic Dual Head Stethoscope

Liberty



Lightweight aluminium stethoscope combines outstanding performance with distinctive styling.

A very economical and well designed unit.

The latex-free dual head stethoscope features a lightweight polished aluminium dual head chest piece and flexible 56cm PVC tubing. Overall length is 76cm. The chrome plated brass binaurals can rotate for a comfortable fit, and come complete with PVC ear tips.

Features:

- Light weight aluminium reversible head chest piece
- Adjustable chrome plated binaural
- Flexible 56cm PVC tubing
- Non chill ring for patient comfort
- Boxed
- Warranty: 12 Month Manufacturer's Warranty
- Stethoscope Length: 76cm



Product Code	Description	
BSCBK	DUAL HEAD STETHOSCOPE BLACK	
BSCR	DUAL HEAD STETHOSCOPE RED	
BSCT	DUAL HEAD STETHOSCOPE TEAL	



Prepared by: Paragon Medical Ltd trading as LR Instruments IRD 113 338 903
 PO Box 145 Greenhithe Auckland NZ 0632
 Contact Details: P 0800 150 030 | F 0800 150 031 | E info@lrinstruments.com.au
 File Name: BSC_Basic-Dual-Head-Steth_v2-0

First Prepared: 21/01/15 Eevien Tan
 Release Approval: 21/01/15 Lee Carter
 Version Control: V2.0 Page 1 | 1