

# Pinard Fetal Stethoscope

## Liberty



Your essential tool for midwifery practice!

The Pinard stethoscope - a trumpet-shaped instrument which can be placed on the maternal abdomen over the fetal chest to hear the fetal heart sounds. Also called a fetal or monoaural stethoscope.

**Features:**

- Cone shaped for exceptional sound response
- BSF Plastic Length: 125mm
- BSFA Aluminium Length: 145mm
- Both Width: 50mm
- No sharp edges
- Warranty: 12 Month Manufacturer's Warranty



BSF



BSFA

Code	Description
BSF	PLASTIC FETAL STETHOSCOPE - ECONOMY
BSFA	ALUMINIUM FETAL STETHOSCOPE - ECONOMY



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107  
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134  
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au  
 File Name: BSF\_Fetal-Steth\_v2-0

First Prepared: 21/01/15 Eevien Tan  
 Release Approval: 21/01/15 Lee Carter  
 Version Control: V2.0 Page 1 | 1