

# Basic Hand-Held Sphygmomanometer

## One-Handed Aneroid

### Features

- Palm type sphygmomanometer
- Latex free
- Designed for convenience in both left-handed and right-handed use
- Large 55cm dial
- CE0123 Approved
- Includes, Adult Cuff, Adult Latex-Free Bladder, Bulb with Intake Valve and a Nylon Carry Case



Complete with a removeable protective ring and nylon zippered carry case



BSOHAV  
Violet



BSOHARB  
Royal Blue

Product Code	Description
BSOHARB	BASIC SPHYG ANEROID ONE-HAND ROYAL BLUE
BSOHAS	BASIC SPHYG ANEROID ONE-HAND SILVER
BSOHAV	BASIC SPHYG ANEROID ONE-HAND VIOLET

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107  
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134  
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au  
File Name: BSOHA\_Sphyg-One-Hand-Aneroid\_v2-1



First Prepared: 14/01/15 Eevien Tan  
Release Approval: 14/01/15 Lee Carter  
Version Control: V2.1 Page 1 | 1