



WAVY FLEXITEM FOR USE INSIDE LONG INSTRUMENTS



A new Wavy Flexistem developed for cleaning the inside of very long instruments where pushing would normally cause buckling.

The wavy design increases its strength enabling its use in very long equipment.

Comes boxed in a 5 metre coil which can be cut to size. Either pull to length and cut or remove from box and place in the Flexistem dispenser FLEX025.



Code
CFLEX007A

Description
FLEXITEM WAVY BRISTLE STEM 5 METRE COIL

All images are included for reference only.
Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes.

E&OE First Published: 1 June 2009
Last Updated: 7 November 2011
File Name: info_CFLEX007A_FLEXITEM-WAVY

