



FLEXISTEM MICROCLENE

FLEXIBLE CLEANING OF SMALL BORE ITEMS



A development of Flexistem, Microclene is more brush like and finer. FLEX021 will clean down to very small diameters and is ideal for Zoellner suction tubes.

Microclene is a fine bore cleaning system for instruments with a very small bore. It can be used to safely dislodge deposits during the washing process.

Flexistem Microclene is used as a brush ensuring that the final rinse process is very effective. As with all Flexistem products, Microclene can be doubled over to increase the cleaning diameter or to increase the rigidity of the wire. Supplied in lengths of 7.5 metre rolls which can be cut to length with the cutters provided, Microclene may be safely disposed of after use.



CFLEX024

CFLEX022

CFLEX023

CFLEX021



Code	Description
CFLEX021	FLEXISTEM MICROCLENE 3MM SUPER FINE
CFLEX022	FLEXISTEM MICROCLENE 6MM MEDIUM FINE
CFLEX023	FLEXISTEM MICROCLENE 3MM SUPER BRUSH
CFLEX024	FLEXISTEM MICROCLENE 6MM MEDIUM BRUSH
CFLEX025	COMPLETE SET INCLUDING DISPENSER, CUTTER & ONE PACK OF EACH ITEM.
CFLEX026	REFILL SET (ONE PACK OF EACH ITEM)

All images are included for reference only.

Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes.

E&OE First Published: 1 June 2009
 Last Updated: 7 November 2011
 File Name: info_CFLEX021-026_FLEXISTEM-MICROCLENE

