

# CLINIPAK Millenium Blue Tray Liners



## Lightweight easy to use viscose base

- Lightweight and very easy to use
- 100% Viscose based dry laid product
- A similar absorbency to super soft soaker sheets but made of a lighter material
- Matt blue finish to minimise glare
- Ideal for medical use

This 40g/m<sup>2</sup> viscose based dry laid liner is ideal for smaller trays and where economy is important. The absorbent qualities assist in the spread of condensate during autoclaving and help achieve dry packs during processing.

Can be supplied in almost any size!

Being very light in weight it will normally be used inside the tray as illustrated.

Conforms extremely well to tray contours - easy to handle

Please Note: The Clinipak Millenium Blue Tray Liners are not suitable for use in a Sterrad Hydrogen Peroxide Steriliser. This product is a Single-Use item.



Code	Description
CMILB002	MILLENNIUM BLUE TRAY LINER 310 X 410 MM (PK 1000)
CMILB004	MILLENNIUM BLUE TRAY LINER 240 X 250 MM (PK 2400)
CMILB005	MILLENNIUM BLUE TRAY LINER 250 X 500 MM (PK 1200)
CMILB016	MILLENNIUM BLUE TRAY LINER 450 X 610 MM (PK 1000)

### Also available from Clinipak

CMILB012	MILLENNIUM BLUE TRAY LINER 180 X 280 MM (PK 2000)
CMILB019	MILLENNIUM BLUE TRAY LINER 305 X 305 MM (PK 2000)
CMILB020	MILLENNIUM BLUE TRAY LINER 500 X 500 MM (PK 1000)

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107  
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134

Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au

File Name: info\_CMILB\_Millenium-Blue-Tray-Liners



Last Modified: 12/09/2014 LC  
First Prepared: 01/06/2009 LC  
Release Approval: 12/09/2014 KS  
Version Control: 2.0 Page 1 | 1