




CLINIPAK POINT PROTECTORS



Point / Needle Protectors


For protecting chisels and other sharp instruments during sterilisation and also transport and storage. Simple and safe to use. Effective, economical and safe to discard.

Available in 2 widths—30mm & 15mm width slits. 'Scissor' type instruments may be inserted opened and locked open by folding the top cut upwards—ensuring good steam circulation. Totally disposable.



Transparent Point Protectors

For the protection of tips and points where visual identification is helpful. Fully autoclavable and suitable for Steam, Plasma and Ethylene Oxide sterilisation. Instrument can be held open to allow steam penetration. The unique design minimises instrument movement.



Code

CPPT001
CPPT002
CPPT003
CPPT004
CPPT006
CPPT007C
CPPT008C

Description

CLINIPAK POINT PROTECTOR 15MM WIDE (BX 2000)
CLINIPAK POINT PROTECTOR 15MM WIDE (BX 1000)
CLINIPAK POINT PROTECTOR 30MM WIDE (BX 2000)
CLINIPAK POINT PROTECTOR 30MM WIDE (BX 1000)
CLINIPAK NEEDLE PROTECTOR (BX 1000)
CLINIPAK POINT PROTECTOR TRANSPARENT 5MM WIDE (BX 200)
CLINIPAK POINT PROTECTOR TRANSPARENT 30MM WIDE (BX 200)

All images are included for reference only.

Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes.

E&OE First Published: 1 June 2009
Last Updated: 7 November 2011
File Name: info_CPPT_CLINIPAK_POINT-PROTECTORS

