



## CLINIPAK TRAY TAGS, TIES & SPLIT RINGS



A range of specially formulated polymer tags in a range of 16 colours with a large tag available in white. These tags may be written upon with non toxic markers or labelled with bar codes as shown.

Clinipak Tray Tags are manufactured from a special mixture of polymers which provide strength and rigidity and a secure surface for affixing labels or writing upon. Labels with permanent adhesives should remain in position for up to 3 years. The tags autoclave well and will not shed or bleach colour.

Three fixing holes enable almost any fixing combination which may be accomplished with the 2 sizes of ties available.

### Identification and traceability

If Tray Tags are used loose within a tray, x-ray identification may be necessary and this is achieved by fixing stainless steel rings to one corner.

The tags are produced in two sizes with the larger size in white only. Colour coding in this vase can be achieved with colour labels.



*A range of specially formulated polymer tags in a range of 10 colours with a larger tag available in white.*

Code	Description
CTAG001	TRAY TAG 90 X 55 MM WHITE (PK 100)
CTAG011	TRAY TAG 90 X 38 MM WHITE (PK 100)
CTAG012	TRAY TAG 90 X 38 MM TAN (PK 100)
CTAG013	TRAY TAG 90 X 38 MM YELLOW (PK 100)
CTAG014	TRAY TAG 90 X 38 MM ORANGE (PK 100)
CTAG015	TRAY TAG 90 X 38 MM RED (PK 100)
CTAG016	TRAY TAG 90 X 38 MM BLUE (PK 100)
CTAG017	TRAY TAG 90 X 38 MM GREEN (PK 100)
CTAG018	TRAY TAG 90 X 38 MM GREY (PK 100)
CTAG019	TRAY TAG 90 X 38 MM PURPLE (PK 100)
CTAG020	TRAY TAG 90 X 38 MM BLACK (PK 100)
CTAG021	TRAY TAG 90 X 38 MM LIME (PK 100)
CTAG022	TRAY TAG 90 X 38 MM PINK (PK 100)
CTAG023	TRAY TAG 90 X 38 MM TURQUOISE (PK 100)
CTAG024	TRAY TAG 90 X 38 MM MUSTARD (PK 100)
CTAG025	TRAY TAG 90 X 38 MM TERRACOTTA (PK 100)
CTAG026	TRAY TAG 90 X 38 MM ROYAL BLUE (PK 100)



*The tags may be written upon with non-toxic markers or labelled with bar codes as shown.*

Above also available in quantities of 20. To order these packs of 20, please add "/20" to the end of the code (ie. CTAG026/20)

All images are included for reference only.

Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes.

E&OE First Published: 1 June 2009  
Last Updated: 7 November 2011  
File Name: info\_CTAG\_CLINIPAK\_TRAY-TAGS





## CLINIPAK TRAY TAGS, TIES & SPLIT RINGS

### Material

Tray Tags provide easy identification of sets or instruments and available in many colours. Surface is textured to allow the use of a marker pen or adherence of self adhesive labels with bar codes. Three holes provided for fixing to split rings or plastic security ties.

- Manufactured from Glass Filled Reinforced Polyamide (PA6), free from latex or latex derivatives
- Glass Content 30 %
- Density 1.36 g/cm<sup>3</sup>
- Water Absorption 6.5 %
- Tensile Strength 170 N/mm<sup>2</sup>
- Elongation at Break 3 %
- Melt Temp. 220 oC
- Split Rings
- Manufactured in stainless steel at 24 mms diameter
- X-Ray detectable
- Plastic Security Ties
- Manufactured from Nylon 66
- Available in Blue (3.6mm width) and clear (2.5mms width), each 150 mms long
- Tray Tags, split rings and plastic security ties can be sterilised by Steam, Ethylene Oxide and Plasma

*Two thicknesses of ties are available for fixing purposes and the three holes in the tags provide for all fixing systems.*



*Autoclavable plastic security ties are available in packets of 100 in either blue or clear.*



*Additional coding is possible by the use of autoclavable self adhesive colour labels size 65 mm x 25 mm.*



*Optional stainless steel rings are available to provide x-ray identification if the tags are used loose in the tray.*



### Code

CLAB002  
CRIN001  
CTIE001  
CTIE002

### Description

SELF ADHESIVE GREEN LABEL 65 X 25MM (PK 180)  
TRAY TAG SPLIT RING 24MM STAINLESS STEEL (PK 100)  
BLUE PLASTIC TIE 3.6 MM (PK 100)  
CLEAR PLASTIC TIE 2.5 MM (PK 100)

### Discontinued Items

CLAB001 SELF ADHESIVE RED LABEL (DISCONTINUED)  
CLAB003 SELF ADHESIVE YELLOW LABEL (DISCONTINUED)  
CLAB004 SELF ADHESIVE BLUE LABEL (DISCONTINUED)  
CLAB005 SELF ADHESIVE WHITE LABEL (DISCONTINUED)

All images are included for reference only.

Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes.

