

Digital Clinical Thermometers

Digital Thermometer Flexible Tip



Features:

- Flexible, soft probe
- Display range: 32.0 to 42.0°C
- Accuracy: $\pm 0.1^{\circ}\text{C}$
- Minimum scale: 0.1
- Measurement time:
 - 60 \pm 10 seconds oral
 - 100 \pm 20 seconds underarm
- Beeper function
- Auto shut-off
- Capacity: 1.5V button battery (LR/SR-41)
- Battery life 100 hours +
- Memory: last measuring reading
- Size: 126 x 18 x 10mm
- LCD: 20 x 7mm
- Net weight: 12g
- Complete with storage case



| Code | Description |
|-------|---|
| CTDFT | THERMOMETER DIGITAL CLINICAL FLEXIBLE TIP |

| | |
|------------------|--|
| Prepared by: | Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134 |
| Contact Details: | P (03) 9870 7066 F (03) 9870 1638 E info@lrinstruments.com.au |
| File Name: | CTDFT_Thermometer-Flexible-Tip_v2-0 |



| | | |
|-------------------|----------|------------|
| First Prepared: | 13/11/14 | Eevien Tan |
| Release Approval: | 17/11/14 | Lee Carter |
| Version Control: | V2.0 | Page 1 1 |