

Digital Clinical Thermometers

Digital Thermometer Rapid Read



Features:

- Display range: 32.0 - 42.0° C
- Accuracy: 0.1° C
- Min. Scale: 0.1
- Measurement time: 10 seconds average
- Beeper function when completed
- Auto shut-off
- Battery: 1.5V button battery
- Size: 134 x 22 x 14mm
- LCD: 20 x 7mm
- Net Weight: 12.5g
- Memory: Last measuring reading
- Complete with storage case



Code	Description
CTDR	THERMOMETER DIGITAL RAPID READ

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: CTDR_Thermometer-Rapid-Read_v2-0



First Prepared: 13/11/14 Eevien Tan
 Release Approval: 17/11/14 Lee Carter
 Version Control: V2.0 Page 1 | 1