

Digital Clinical Thermometers

Digital Thermometer Rigid Tip



Features:

- Display range: 32.0 to 42.0°C
- Accuracy: $\pm 0.1^{\circ}\text{C}$
- Minimum scale: 0.1
- Measurement time:
 - 60 \pm 10 seconds oral
 - 100 \pm 20 seconds underarm
- Beeper function
- Auto shut-off
- Battery: 1.5V button battery (LR/SR-41)
- Battery life 100 hours +
- Memory: last measuring reading
- Size: 126 x 18 x 10mm
- Net weight: 10g
- Complete with storage case



Code	Description
CTD	DIGITAL THERMOMETER RIGID TIP

Prepared by:	Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details:	P (03) 9870 7066 F (03) 9870 1638 E info@lrinstruments.com.au
File Name:	CTD_Thermometer-Rigid-Tip_v2-0



First Prepared:	13/11/14	Eevien Tan
Release Approval:	17/11/14	Lee Carter
Version Control:	V2.0	Page 1 1