

# Liberty Range of Thermometers



## Digital Thermometer Rigid Tip

### Features:

- Display range: 32.0 to 42.0°C
- Accuracy:  $\pm 0.1^{\circ}\text{C}$
- Minimum scale: 0.1
- Measurement time:
  - 60  $\pm$  10 seconds oral
  - 100  $\pm$  20 seconds underarm
- Beeper function
- Auto shut-off
- Battery: 1.5V button battery (LR/SR-41)
- Battery life 100 hours +
- Memory: last measuring reading
- Size: 126 x 18 x 10mm
- Net weight: 10g
- Complete with storage case

Ordering Code: CTD  
Digital Thermometer Rigid Tip



Code	Description
CTD	DIGITAL THERMOMETER RIGID TIP

## Digital Thermometer Flexible Tip

### Features:

- Flexible, soft probe
- Display range: 32.0 to 42.0°C
- Accuracy:  $\pm 0.1^{\circ}\text{C}$
- Minimum scale: 0.1
- Measurement time:
  - 60  $\pm$  10 seconds oral
  - 100  $\pm$  20 seconds underarm
- Beeper function
- Auto shut-off
- Capacity: 1.5V button battery (LR/SR-41)
- Battery life 100 hours +
- Memory: last measuring reading
- Size: 126 x 18 x 10mm
- LCD: 20 x 7mm
- Net weight: 12g
- Complete with storage case

Ordering Code: CTDF  
Digital Thermometer Flexible Tip



Code	Description
CTDF	THERMOMETER DIGITAL CLINICAL FLEXIBLE TIP



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107  
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134  
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au  
File Name: CT\_Thermometers-Liberty\_v2-0

First Prepared: 13/11/14 Eevien Tan  
Release Approval: 17/11/14 Lee Carter  
Version Control: V2.0 Page 1 | 2

# Liberty Range of Thermometers



## Digital Thermometer Rapid Read

### Features:

- Display range: 32.0 - 42.0° C
- Accuracy: 0.1° C
- Min. Scale: 0.1
- Measurement time: 10 seconds average
- Beeper function when completed
- Auto shut-off
- Battery: 1.5V button battery
- Size: 134 x 22 x 14mm
- LCD: 20 x 7mm
- Net Weight: 12.5g
- Memory: Last measuring reading
- Complete with storage case

Ordering Code: CTDR  
Digital Thermometer Rapid Read



Code	Description
CTDR	THERMOMETER DIGITAL RAPID READ

## Tympanic Ear Thermometer

### Features:

- Safe and Hygienic
- Easy-to-read, Three-digit LCD display
- Auto shut-off
- Memory: 10 recording
- Measuring time: 1 second with beeper function
- Battery: 2 x AAA

Ordering Code: CTT  
Tympanic Ear Thermometer



Code	Heading
CTT	TYMPANIC INFARED EAR THERMOMETER
CTTC20	DISPOSABLE COVERS (BOX 20)

## Stubby Bulb Clinical Thermometer

### Features:

- Accurate
- Easy to read
- Mercurial
- Complete with storage case

Ordering Code: CTSBC  
Stubby Bulb Thermometer



Code	Heading
CTSBC	STUBBY BULB CLINICAL THERMOMETER



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107  
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134  
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au  
File Name: CT\_Thermometers-Liberty\_v2-0

First Prepared: 13/11/14 Eevien Tan  
Release Approval: 17/11/14 Lee Carter  
Version Control: V2.0 Page 2 | 2