

Basic Dissecting Kit

A quality dissection kit, geared specifically for school, college students and teachers, are a great learning and exploring tool. Containing instruments necessary to perform complete dissection in a vinyl zippered wallet. This dissecting kit is designed to meet the requirements of students in middle/high school right through to university level anatomy courses.

A dissection lab will require use of dissecting kits. We offer a range of dissecting kits that are tailored to meet specific lab requirements. These dissecting kits will generally be geared towards either plant or animal dissection.

Each academic level requires different dissection kit. Smaller dissecting kits (basic dissection kit) will be for elementary age, which will focus on small specimens such as flowers or earthworms. Larger, more advanced dissecting kits will be for dissecting frog, fetal pigs and larger preserved animals.

If you think you need more tools you can order from our large collection of dissecting tools along with a dissecting kit(s) that is more suitable to your needs.

Complete with

- 1 Seeker
- 2 Mounted Dissecting Needles
- 1 Scalpel Blade Handle
- 1 Blunt Forceps
- 1 Pointed Forceps
- 1 pair 12.5cm SH/BL Surgical Scissors
- 3 Scalpel Blades
- 1 pair 10cm Iris Dissecting Scissors
- Supplied in a zippered wallet



Code	Description
DKBASIC	BASIC DISSECTING KIT

Prepared by:	Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details:	P (03) 9870 7066 F (03) 9870 1638 E info@lrinstruments.com.au
File Name:	DKBASIC_Dissecting-Kit_v4-0



First Prepared:	14/11/14	Eevien Tan
Release Approval:	14/11/14	Lee Carter
Version Control:	V4.0	Page 1 1

Disclaimer: All images are included for reference only. Features and specifications are subject to change without prior notice.
Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes. E & OE
Visit our website at: <http://www.lrinstruments.com.au>