

UltraTec Pocket Doppler

PD1CV

Vascular Pocket Doppler

The UltraTec PD1CV is a small, lightweight, high sensitivity pocket Doppler with fixed vascular probes. It is perfectly suited as an entry level doppler for the General practitioner or emergency nurse as well as Podiatrists, Chiropodists or vascular specialists for pressure index or routine vascular detection. Choosing a fixed 5 or 8 MHz vascular probe provides an ideal tool for use by a General Practitioner for routine vascular investigation, in ER for vessel location or for monitoring pulses and systolic blood pressure and also by the specialist in diabetes or dermatology.

Audio signals are available to the user via the built-in loudspeaker or optional headset with output level controlled by rotary volume control.

The unit is powered by a standard 9V battery which provides over 9 hours of use. The PD1CV meets all current EEC standards for medical equipment.

Features

- Lightweight unit and Probe
- Easy to clean
- 5 or 8MHz vascular probes
- Battery saver system
- 3 year warranty



Code	Description
DPD1CV5	PD1CV VASCULAR DOPPLER WITH 5 MHZ PROBE
DPD1CV8	PD1CV VASCULAR DOPPLER WITH 8 MHZ PROBE
DPD1V5	PD1V VASCULAR DOPPLER WITH 5 MHZ PENCIL PROBE
DPD1V8	PD1V VASCULAR DOPPLER WITH 8 MHZ PENCIL PROBE

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: DPD1CV_Doppler-Ultratec_Vascular_v1-2



First Prepared: 14/11/14 Eevien Tan
 Release Approval: 14/11/14 Lee Carter
 Version Control: V1.2 Page 1 | 1