

UltraTec Pocket Doppler PD1

Pocket Fetal Doppler

The UltraTec PD1 is a small, lightweight, high sensitivity pocket Doppler. Its 2MHz fetal Doppler probe is ideal for use by General Practitioner, Obstetrician or Midwife for routine antenatal fetal detection. The fetal signal is available to the user via the built in loudspeaker or optional headset with output level controlled by edge mounted volume control. The fetal signals are easily located using the lightweight narrow beam transducer while internal signal processing minimises the noise artefacts.

The unit is powered by a standard 9V battery which provides over 9 hours of use. The PD1 meets all current EEC standards for medical equipment.

Features

- Lightweight ABS unit and probes (295g typically including transducer)
- Easy to clean ABS housing
- 2MHz fixed transducer
- Battery saver system with 9 hour life
- Audio output for optional headphones
- Includes soft carry bag
- 3 year warranty



Code	Description
DPD1	PD1 FETAL DOPPLER 2 MHZ PROBE
DPDAH	DOPPLER AUDIO HEADSET
DPD1SC	DOPPLER SPARE SOFT CASE

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: DPD1_Doppler-Ultratec-Fetal_v1-2



First Prepared: 14/11/14 Eevien Tan
 Release Approval: 14/11/14 Lee Carter
 Version Control: V1.2 Page 1 | 1