

Liberty Penlight Torches



Economy Penlight Torch

Features

- Central push button on/off switch
- Focusing lens
- Pocket clip
- Sealed concentric beam bulb
- Takes 2 x AAA size batteries (not included)
- Standard illumination 2.5V
- Enamel coated aluminium case


CODE	DESCRIPTION
DTEB	ECONOMY PENLIGHT TORCH BLACK
DTEBCS	ECONOMY PENLIGHT TORCH BLACK IN CLAM SHELL
DTEW	ECONOMY PENLIGHT TORCH WHITE
DTEWCS	ECONOMY PENLIGHT TORCH WHITE IN CLAM SHELL
DTEPI-MF	ECONOMY PENLIGHT TORCH MCGRATH PINK
DTEPICS-MF	ECONOMY PENLIGHT TORCH MCGRATH PINK IN CLAM SHELL



DTEB



DTEPI-MF




DTEW

Disposable Penlight Torches

The modern way to deliver self focusing, high intensity bright light without the inconvenience of changing batteries and bulbs. Compact and lightweight the Liberty Disposable Penlight Torch provides a clear, crisp light right where it's needed. Long-lasting, disposable, featuring convenient clip-activated switch. Constant use time of 30 minutes.

Features

- Sealed unit with 2 x AAA size batteries included
- Concentric beam bulb with pocket clip/switch
- Completely disposable
- White, sold in six pack



CODE	DESCRIPTION
TDL	LIBERTY DISPOSABLE PENLIGHT TORCH (PK 6)

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: DTE_Liberty_Penlight-Torches_v4-1



First Prepared: 14/11/14 Eevien Tan
 Release Approval: 17/11/14 Lee Carter
 Version Control: V4.1 Page 1 | 2


Liberty Penlight Torches



Luxury Penlight Torch

Features

- Modern design penlight
- Stainless steel casing DTSS
- Enamel coated aluminium casing (DTL)
- Protected recessed bulb on DTL
- Central push button on/off switch
- Sealed concentric beam bulb
- Pocket clip
- Standard illumination 2.5V
- Takes 2 x AAA size batteries (included)
- Weight: 14.3g; Length: 136mm; Width: 12mm



DTLRB



DTLPU



DTLW



DTLB



DTSS

CODE	DESCRIPTION
DTLB	LIBERTY LUXURY PENLIGHT TORCH BLACK
DTLRB	LIBERTY LUXURY PENLIGHT TORCH BLUE
DTLPU	LIBERTY LUXURY PENLIGHT TORCH PURPLE
DTLW	LIBERTY LUXURY PENLIGHT TORCH WHITE
DTSS	LUXURY PENLIGHT TORCH STAINLESS STEEL

Supplied in Clam Shells

DTLBSCS	LIBERTY LUXURY PENLIGHT TORCH BLACK IN CLAM SHELL
DTLRBCS	LIBERTY LUXURY PENLIGHT TORCH BLUE IN CLAM SHELL
DTLPUCS	LIBERTY LUXURY PENLIGHT TORCH PURPLE IN CLAM SHELL
DTLWCS	LIBERTY LUXURY PENLIGHT TORCH WHITE IN CLAM SHELL
DTSSCS	LUXURY PENLIGHT TORCH STAINLESS STEEL IN CLAM SHELL

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: DTE_Liberty_Penlight-Torches_v4-1



First Prepared: 13/11/14 Eevien Tan
 Release Approval: 17/11/14 Lee Carter
 Version Control: V4.1 Page 2 | 2