


Luxury Penlight Torch

Liberty



Features

- Modern design penlight
- Stainless steel casing DTSS
- Enamel coated aluminium casing (DTL)
- Protected recessed bulb on DTL
- Central push button on/off switch
- Sealed concentric beam bulb
- Pocket clip
- Standard illumination 2.5V
- Takes 2 x AAA size batteries (included)
- Weight: 14.3g; Length: 136mm; Width: 12mm



Code	Description
DTLB	LIBERTY LUXURY PENLIGHT TORCH BLACK
DTLRB	LIBERTY LUXURY PENLIGHT TORCH BLUE
DTLPU	LIBERTY LUXURY PENLIGHT TORCH PURPLE
DTLW	LIBERTY LUXURY PENLIGHT TORCH WHITE
DTLS	LIBERTY LUXURY PENLIGHT TORCH SILVER

Supplied in Clam Shells	
DTLBCS	LIBERTY LUXURY PENLIGHT TORCH BLACK IN CLAM SHELL
DTLRBCS	LIBERTY LUXURY PENLIGHT TORCH BLUE IN CLAM SHELL
DTLPUCS	LIBERTY LUXURY PENLIGHT TORCH PURPLE IN CLAM SHELL
DTLWCS	LIBERTY LUXURY PENLIGHT TORCH WHITE IN CLAM SHELL
DTLSCS	LIBERTY LUXURY PENLIGHT TORCH SILVER IN CLAM SHELL

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: DT_Luxury-Penlight-Torches_v2-1



First Prepared: 30/03/14 Eevien Tan
 Release Approval: 30/03/14 Lee Carter
 Version Control: V2.1 Page 1 | 1