

ERKA Pressure Infusion

Pressure infusion devices are an essential item in any hospital. Pressure cuffs in combination with a manometer and inflation system provide for optimal pressure infusions. Suitable for solutions or blood bags in the following sizes: 500 ml and 1000 ml.

Features

- Sturdy nylon cuff in blue, with suspension loop for infusion stands
- Washable up to 60°C
- Precise reading of the level in the bag thanks to a transparent net on the cuff
- Manometer with 39 mm diameter aluminium scale up to 750 mm Hg



Code	Description
EP821 000 00	ERKA PRESSURE INFUSER 500 ML
EP831 000 00	ERKA PRESSURE INFUSER 1000 ML

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: EP8_ERKA-Pressure-Infuser_v2-0



First Prepared: 14/11/14 Eevien Tan
 Release Approval: 17/11/14 Lee Carter
 Version Control: V2.0 Page 1 | 1