

Adult Cuff & Bladder Sets



Features

- Latex Free inflation system
- Index and range and artery indicator label markings
- Adult Range: 29-40cm
- Velcro fastener
- Extra long single tube 80cm in length
- Comes as standard with a nylon adult cuff and latex-free bladder 1-tube



Product Code	Description	
ESCBAKLF	LIBERTY CUFF & BLADDER SET ADULT 1-TUBE LATEX FREE BLACK	
ESCBAPILF-MF	LIBERTY CUFF & BLADDER SET ADULT 1-TUBE LATEX FREE MCGRATH PINK	
ESCBARBLF	LIBERTY CUFF & BLADDER SET ADULT 1-TUBE LATEX FREE ROYAL BLUE	
ESCBVLF	LIBERTY CUFF & BLADDER SET ADULT 1-TUBE LATEX FREE VIOLET	



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: ESCBA_Cuff-Bladder-Set_v2-1

First Prepared: 15/01/15 Eevien Tan
 Release Approval: 15/01/15 Lee Carter
 Version Control: V2.1 Page 1 | 1