

ERKA Vario

Aneroid Sphyg Models



ERKA Vario Desk Model

Features

- High quality precision manometer
- Diameter of scale 150 mm, with mirrored exterior ring
- Scale inclined by 12°
- Housing made of high-quality ABS in signal white/light-grey
- Non-reflective manometer-glass
- Cuff-basket at the back in powder-coated steel
- Inflation bulb with sensitive adjustable screw reduction valve
- Measuring Range: 0 - 300 mm Hg
- Dimensions: W175mm x H175mm x D130mm



Sphyg Erka vario desk model aneroid sphyg

Product Code	Description
EV285DM	ERKA VARIO DESK ANEROID SPHYG WITH ADULT CUFF

ERKA Vario Wall Model

The Vario basic model combines with simple accessories for universal use:

- Wall model with swivel motion
- Rail-mounted model with swivel motion
- Stand/trolley model on 5 castors

Features

- High quality precision manometer
- Diameter of scale 150 mm, with mirrored exterior ring
- Scale inclined by 12°
- Housing made of high-quality ABS in signal white/light-grey
- Non-reflective manometer-glass
- Cuff-basket at the front in powder-coated steel
- Inflation bulb with sensitive adjustable screw reduction valve
- Measuring Range: 0 - 300 mm Hg
- Dimensions: W175mm x H250mm x D40mm



Erka vario wall mounted aneroid sphyg

Product Code	Description
EV288WA	ERKA VARIO WALL ANEROID SPHYG WITH ADULT CUFF

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: EV28_Erka-Vario-Aneroid-Sphyg_v3-1



First Prepared: 15/01/15 Eevien Tan
 Release Approval: 15/01/15 Lee Carter
 Version Control: V3.1 Page 1 | 2

ERKA Vario

Aneroid Sphyg Models



ERKA Vario Trolley Model

Features

- High quality precision manometer
- Diameter of scale 150mm, with mirrored exterior ring
- Scale inclined by 12°
- Housing made of high quality ABS in signal white/light-grey
- Non-reflective manometer glass
- Separate cuff basket in powder-coated steel
- Inflation bulb with sensitive adjustable screw reduction valve



Erka vario trolley/mobile
Aneroid sphyg

Product Code	Description
EV288ML	ERKA VARIO MOBILE ANEROID SPHYG WITH LIGHT WEIGHT BASE
EV288MM	ERKA VARIO MOBILE ANEROID SPHYG WITH MEDIUM WEIGHT BASE
EV288MH	ERKA VARIO MOBILE ANEROID SPHYG WITH HEAVY WEIGHT BASE

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: EV28_Erka-Vario-Aneroid-Sphyg_v3-1



First Prepared: 15/01/15 Eevien Tan
Release Approval: 15/01/15 Lee Carter
Version Control: V3.1 Page 2 | 2