

Safety Glasses

Liberty



Excellent optical quality, these modern, stylish Liberty Safety Glasses are ideal for doctors and all medical personnel. They offer great eye protection with a wide field of vision with built-in side shields making a comfortable face hugging fit.

Features

- Ideal for doctors and all medical personnel
- Frameless safety spectacles for unrestricted vision in all directions
- Soft non-allergenic Anti-slip temple inserts
- Soft Nose-bridge for maximum comfort
- Built-in side shields for added protection
- Low energy impact, UV Filtering, Sun glare Protection
- Tested and Approved to CE EN166, 170, 172 & ANSI Z87.1 standards



Safety glasses or goggles are a form of protective eye wear that usually enclose or protect the eye area in order to prevent particulates, water or chemicals from striking the eyes. Safety glasses should be worn to prevent flying particles from damaging the eyes.

Safety glasses are usually made with shatter-resistant plastic lenses to protect the eye from flying debris. Although safety lenses may be constructed from a variety of materials of various impact resistance, certain standards suggest that they maintain a minimum 1 millimeter thickness at the thinnest point, regardless of material. Safety glasses can vary in the level of protection they provide. For example, those used in medicine may be expected to protect against blood splatter while safety glasses in a factory might have stronger lenses and a stronger frame with additional shields at the temples.

CODE	DESCRIPTION
GSL	LIBERTY SAFETY GOGGLES

Prepared by:	Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details:	P (03) 9870 7066 F (03) 9870 1638 E info@lrinstruments.com.au
File Name:	GSL_Safety-Glasses_v2-0



First Prepared:	14/11/14	Eevien Tan
Release Approval:	17/11/14	Lee Carter
Version Control:	V2.0	Page 1 1