

Safety Glasses

Safety glasses or goggles are a form of protective eye wear that usually enclose or protect the eye area in order to prevent particulates, water or chemicals from striking the eyes. Safety glasses should be worn to prevent flying particles from damaging the eyes.

With excellent optical quality, these quality safety glasses are ideal for doctors and all medical personnel. They offer great eye protection with the side shields and a wide field of vision and offer UV protection to 400 nanometers. With an adjustable temple arm offering positive locks you can achieve a comfortable face hugging fit.

Features

- Adjustable arms
- Ideal for doctors and all medical personnel
- Eye protection
- Wide field of vision
- Excellent optical quality
- U.V. protection to 400 nanometers
- Side shields
- Complies to AS337 - 1992 (low impact)



CODE	DESCRIPTION
GSS	SAFETY GOGGLES

Prepared by:	Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details:	P (03) 9870 7066 F (03) 9870 1638 E info@lrinstruments.com.au
File Name:	GSS_Safety-Glasses_v2-0



First Prepared:	14/11/14	Eevien Tan
Release Approval:	17/11/14	Lee Carter
Version Control:	V2.0	Page 1 1