

Hot & Cold Packs

Instant Cold Compress - UREA

The Liberty Instant Cold Compress temporarily relieves minor pain and swelling for sprains, aches and sore joints. It is conveniently disposable with no pre-chilling required for quick, effective relief. Product becomes cold on its own after activation.

The Liberty Instant Cold Compress is manufactured with urea instead of ammonium nitrate. Not only is this type of instant ice safer, but it is better for the environment. Convenient and portable, for use whenever and wherever cold therapy is needed. Store in a cool, dry place.

Usage

Simply squeeze the pouch. The contents become immediately cold and remain cold for 14-40 minutes and achieves temperatures as low as -2°C

Details

- Size: H228mm x W152mm
- Weight: 260g
- Pack Size: 40 per carton



Code	Description
INTICEUREA	INSTANT COLD COMPRESS (BOX 40)

Prepared by:	Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details:	P (03) 9870 7066 F (03) 9870 1638 E info@lrinstruments.com.au
File Name:	INTCC_Cold-Compress_v3-0



First Prepared:	18/11/14	Eevien Tan
Release Approval:	18/11/14	Lee Carter
Version Control:	V3.0	Page 1 1