



LAWTON INSTRUMENT CLEANING BRUSHES



Lawton instrument cleaning brushes with stainless steel handles.



Please Note:

The Lawton Cleaning Brushes LA07-2030 and LA07-2032 are made out of Nylon and for this reason they are not autoclavable. The other Lawton branded brushes have natural bristles and can be steam sterilised at 121°C maximum. As all the brushes could be heavily contaminated after use, Lawton recommend that they be purchased/sold as a disposable item.

Code	Description
LA07-2038	CLEANING BRUSH 330 MM (BRUSH LENGTH 130 X BRUSH DIAMETER 50 MM)
LA07-2028	CLEANING BRUSH 215 MM (BRUSH LENGTH 100 X BRUSH DIAMETER 28 MM)
LA07-2026	CLEANING BRUSH 215 MM (BRUSH LENGTH 95 X BRUSH DIAMETER 25 MM)
LA07-2024	CLEANING BRUSH 215 MM (BRUSH LENGTH 95 X BRUSH DIAMETER 20 MM)
LA07-2022	CLEANING BRUSH 215 MM (BRUSH LENGTH 65 X BRUSH DIAMETER 15 MM)
LA07-2020	CLEANING BRUSH 215 MM (BRUSH LENGTH 60 X BRUSH DIAMETER 13 MM)
LA07-2036	CLEANING BRUSH 265 MM (BRUSH LENGTH 65 X BRUSH DIAMETER 8 MM)
LA07-2034	CLEANING BRUSH 285 MM (BRUSH LENGTH 75 X BRUSH DIAMETER 7 MM)
LA07-2030	NYLON CLEANING BRUSH 265 MM (BRUSH LENGTH 65 X BRUSH DIAMETER 3 MM)
LA07-2032	NYLON CLEANING BRUSH 430 MM (BRUSH LENGTH 65 X BRUSH DIAMETER 3 MM)

All images are included for reference only.
Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes.

E&OE First Published: 1 June 2009
Last Updated: 10 November 2011
File Name: info_LA07_LAWTON_CLEANING-BRUSHES

