

Laryngoscope

Fibre Optic



Features

- Autoclavable to 134°C
- Sure grip handle
- Takes 2 x 'C' size batteries (not included)
- Satin non-glare finish
- Cold fibre optic light
- Fully compatible with all fibre optic hook-on handles that conform to ASTM (F965)



Code	Description
------	-------------

(NOTE: MCINTOSH BLADE IS CURVED)

LLBFM0	MCINTOSH BLADE FIBRE OPTIC NO. 0
LLBFM1	MCINTOSH BLADE FIBRE OPTIC NO. 1
LLBFM2	MCINTOSH BLADE FIBRE OPTIC NO. 2
LLBFM3	MCINTOSH BLADE FIBRE OPTIC NO. 3
LLBFM4	MCINTOSH BLADE FIBRE OPTIC NO. 4
LLBFM5	MCINTOSH BLADE FIBRE OPTIC NO. 5
LBSFM3	LARYNGOSCOPE SET WITH 3 BLADES AND HANDLE IN CARRY CASE

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107

Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134

Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au

File Name: LLB_Laryngoscope-Fiber-Optic_v3-2



First Prepared: 12/03/15

Eevien Tan

Release Approval: 12/03/15

Lee Carter

Version Control: V3.2

Page 1 | 2

Disclaimer: All images are included for reference only. Features and specifications are subject to change without prior notice.

Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes. E & OE

Visit our website at: <http://www.lrinstruments.com.au>

Laryngoscope

Fibre Optic



Code	Description
------	-------------

(NOTE: MILLER BLADE IS STRAIGHT)

LLBFMI00	MILLER BLADE FIBRE OPTIC NO. 00 (NOT PICTURED)
LLBFMI0	MILLER BLADE FIBRE OPTIC NO. 0
LLBFMI1	MILLER BLADE FIBRE OPTIC NO. 1
LLBFMI2	MILLER BLADE FIBRE OPTIC NO. 2
LLBFMI3	MILLER BLADE FIBRE OPTIC NO. 3
LLBFMI4	MILLER BLADE FIBRE OPTIC NO. 4
LLFH	FIBRE OPTIC HANDLE
LLFHAA	FIBRE OPTIC HANDLE FOR AA SIZE BATTERIES
LLHFS	FIBRE OPTIC HANDLE STUBBY
LLFB	HALOGEN BULB FOR FIBRE OPTIC HANDLE
LLB	REPLACEMENT OPTICAL BUNDLE
LLHFOLA	LIGHT ASSEMBLY TO SUIT LLFH

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: LLB_Laryngoscope-Fiber-Optic_v3-2



First Prepared:	12/03/15	Eevien Tan
Release Approval:	12/03/15	Lee Carter
Version Control:	V3.2	Page 2 2

Disclaimer: All images are included for reference only. Features and specifications are subject to change without prior notice.
 Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes. E & OE
 Visit our website at: <http://www.lrinstruments.com.au>