

Laryngoscope Set

Fibre Optic



Features

- Stainless steel
- Autoclavable to 134°C
- Cold fibre optic light
- MacIntosh Blades in three (3) sizes: Size 1, Size 2, and Size 3
- Fully compatible with all fibre optic hook-on handles that conform to ASTM (F965)
- Halogen bulb
- Sure grip handle
- Satin non-glare finish
- Takes 2 x 'C' size batteries (not included)
- Carry case



Code	Description
LLSFM3	LARYNGOSCOPE SET FIBRE-OPTIC WITH 3 BLADES & CARRY CASE

Prepared by:	Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details:	P (03) 9870 7066 F (03) 9870 1638 E info@lrinstruments.com.au
File Name:	LLSFM3_Laryngoscope-Set_v4-0



First Prepared:	18/11/14	Eevien Tan
Release Approval:	18/11/14	Lee Carter
Version Control:	V4.0	Page 1 1

Disclaimer: All images are included for reference only. Features and specifications are subject to change without prior notice.
Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes. E & OE
Visit our website at: <http://www.lrinstruments.com.au>