

3.5V Coaxial Ophthalmoscope

Rechargeable

Liberty Diagnostic



Coaxial optics enable easy entry into undilated pupils for more accurate diagnoses, while two filter selections provide optimum light levels for entering healthy eyes and those with cataracts.

Features

- Precise lamp filament and optical alignment
- Rechargeable handle - 120 minutes run time
- 3.5V halogen lamp provides light output for true tissue colour and consistent, long-lasting illumination
- 28 focusing lenses with a range of $\pm 1D$ to $+40D$ diopters
- 18 unique aperture/filter combinations for greater versatility

Set Includes

- 1 x 3.5V Ophthalmoscope Head
- 1 x Rechargeable Handle (120 minutes run time)
- 1 x Desk Charger Single Cell
- 1 x Power Adapter/Transformer



6 Standard Apertures



Cobalt blue filter, Micro spot, small spot, large spot, Slit, Fixation aperture.

With Red-Free Filter



Cobalt blue filter, Micro spot, small spot, large spot, Slit, Fixation aperture.

WARRANTY:

1 Year Manufacturers Warranty excluding bulbs and globes. The Rechargeable Handle (with 120 minutes run time) has a 2 Year Warranty. There is no internal battery. A new handle will need to be purchased once the handle no longer works.

Order Code	Description
LOPCRS	3.5V OPHTHALMOSCOPE RECHARGEABLE SET
LOPCRSG	3.5V OPHTHALMOSCOPE SPARE BULB
LOHR	3.5V RECHARGEABLE HANDLE ONLY
LOPCRHO	3.5V OPHTHALMOSCOPE HEAD ONLY
LODCS	DESK CHARGER SINGLE CELL WITH TRANSFORMER
LOT	POWER ADAPTER / TRANSFORMER ONLY

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: LOPCRS_Ophthalmoscope-Rechargeable_v3-0



First Prepared: 19/11/14 Eevien Tan
Release Approval: 19/11/14 Lee Carter
Version Control: V3.0 Page 1 | 1