

Mini Ophthalmoscope

Liberty Diagnostic



The Liberty Ophthalmoscope fits into your pocket!

It's easy to use, lightweight and perfect for examination of the eye and it's background. Take your Liberty Ophthalmoscope with you everywhere you go throughout the day. Available in black only.

Features

- Halogen light for true tissue colour and consistent, long-lasting illumination
- Apertures include large and small spot, semi-circle and red-free filter
- Rubber brow rest prevents scratching of eye-glasses
- Pocket clip with integral on/off slide switch for ambidextrous use
- Supplied in a hard protective carry case

Technical Data


Correcting Lens Range:

+ 1 | 2 | 3 | 4 | 6 | 8 | 10 | 15 | 20 | 0 Diopter

- 1 | 2 | 3 | 4 | 6 | 8 | 10 | 15 | 20 Diopter

Aperture Wheel: Large spot, Small spot, Hemi spot, Red-free filter
 Bulb: 2.5V 1.6W Halogen lamp
 Batteries: 2 x "AA" batteries
 or 2 x "AA" NiMH Rechargeable batteries
 (batteries not included)

Apertures



Order Code **Description**
 LOPM LIBERTY MINI OPHTHALMOSCOPE

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: LOPM_Ophthalmoscope-Mini_v3-0



First Prepared: 19/11/14 Eevien Tan
 Release Approval: 19/11/14 Lee Carter
 Version Control: V3.0 Page 1 | 1