

Fibre Optic Otoscope

Rechargeable

Liberty Diagnostic



The Liberty Fibre Optic Otoscope 3.5V has a sealed system for pneumatic otoscopy and accurate diagnosis of otitis media and a swivelling wide-angle viewing lens.

Features

- Rechargeable Fibre Optic Otoscope 3.5V
- Swivelling wide-angle magnifying lens
- Reusable ear tips
- Rechargeable handle, 120 minutes run time
- Compatible with Welch Allyn disposable ear tips

The Otoscope is used to view the tympanic membrane and external ear canal to diagnose outer and middle ear pathologies. It is the most frequently used device in ENT.

Set Includes

- 1 x 3.5V Otoscope Head
- 1 x Rechargeable Handle, 120 minutes run time
- 1 set (2.75mm and 4.25mm, 5 pcs each)
Reusable Specula (Ear Tips)
- 1 x Desk Charger Single Cell
- 1 x Power Adapter/Transformer



WARRANTY:
1 Year Manufacturers Warranty excluding bulbs and globes. The Rechargeable Handle (with 120 minutes run time) has a 2 Year Warranty. There is no internal battery. A new handle will need to be purchased once the handle no longer works.

Order Code	Description
LOTCSR	3.5V OTOSCOPE RECHARGEABLE SET
WA52432-U	DISPOSABLE EAR TIPS 2.75MM (PK 850)
WA52432-34	DISPOSABLE EAR TIPS 2.75MM (PK 34)
WA52434-U	DISPOSABLE EAR TIPS 4.25MM (PK 850)
WA52434-34	DISPOSABLE EAR TIPS 4.25MM (PK 34)
LOTCRSG	3.5V OTOSCOPE SPARE BULB
LOTCRHO	3.5V OTOSCOPE HEAD ONLY
LOHR	3.5V RECHARGEABLE HANDLE ONLY
LODCS	DESK CHARGER SINGLE CELL WITH TRANSFORMER
LOT	POWER ADAPTER/TRANSFORMER ONLY

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: LOTCRS_Otoscope-Fibre-Optic-Rechargeable_v3-0



First Prepared: 19/11/14 Eevien Tan
Release Approval: 19/11/14 Lee Carter
Version Control: V3.0 Page 1 | 1