

Fibre Optic Otoscope

Liberty Diagnostic



High quality, yet simply designed fibre optic otoscope.
A classic instrument for the physician who requires fibre optic illumination. Top "Liberty" quality at a reasonable price!

Features

- New improved switch and bayonet fitting for a more secure attachment of specula
- 2.5x Swivel Magnifying Lens
- Reusable Ear Tips
- Battery handle (2 x C batteries)
- Top Liberty quality at a reasonable price
- Compatible with Heine disposable ear tips

The Otoscope is used to view the tympanic membrane and external ear canal to diagnose outer and middle ear pathologies. It is the most frequently used device in ENT.



Order Code	Description
LOT C	OTOSCOPE FIBRE OPTIC C
LOTETD2.5	DISPOSABLE EAR TIPS 2.5MM (PK 50)
LOTETD4	DISPOSABLE EAR TIPS 4MM (PK 50)

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: LOTC_Otoscope-Fibre-Optic_v3-0



First Prepared:	19/11/14	Eevien Tan
Release Approval:	19/11/14	Lee Carter
Version Control:	V3.0	Page 1 1