

# Otoscope/Ophthalmoscope Set

## Rechargeable

### Liberty Diagnostic



#### Otoscope Features

The Liberty Fibre Optic Otoscope 3.5V has a sealed system for pneumatic otoscopy and accurate diagnosis of otitis media and a swivelling wide-angle viewing lens allowing for instrumentation under magnification.

- Rechargeable 3.5V Fibre Optic Otoscope
- 2.0 x Swivelling Wide-Angle Magnifying Lens
- Reusable Ear Tips
- Rechargeable handle, 120 minutes run time
- Compatible with Welch Allyn disposable ear tips

#### Ophthalmoscope features

Coaxial optics enable easy entry into undilated pupils for more accurate diagnoses, while two filter selections provide optimum light levels for entering healthy eyes and those with cataracts.

- Rechargeable 3.5V Coaxial Ophthalmoscope
- Precise lamp filament and optical alignment
- Provides light output for true tissue colour and consistent, long-lasting illumination
- 28 focusing lenses with a range of  $\pm 1D$  to  $+40D$  diopter's
- 18 unique aperture/filter combinations for greater versatility
- Rechargeable handle, 120 minutes run time

#### Set Includes

- 1 x 3.5V Otoscope Head
- 1 x 3.5V Ophthalmoscope Head
- 2 x Rechargeable Handles, 120 minutes run time
- 1 set (2.75mm and 4.25mm, 5 pcs each)  
Reusable Specula (Ear Tips)
- 1 x Desk Charger Double Cell
- 1 x Power Adapter/Transformer



#### WARRANTY:

1 Year Manufacturers Warranty excluding bulbs and globes. The Rechargeable Handle (with 120 minutes run time) has a 2 Year Warranty. There is no internal battery. A new handle will need to be purchased once the handle no longer works.

Order Code	Description
LOTPCSET	3.5V OTOSCOPE/OPHTHALMOSCOPE RECHARGEABLE SET
WA52432-U	DISPOSABLE EAR TIPS 2.75MM (PK 850)
WA52432-34	DISPOSABLE EAR TIPS 2.75MM (PK 34)
WA52434-U	DISPOSABLE EAR TIPS 4.25MM (PK 850)
WA52434-34	DISPOSABLE EAR TIPS 4.25MM (PK 34)
LOTCRSG	SPARE GLOBE FOR 3.5V OTOSCOPE
LOPCRSG	SPARE GLOBE FOR 3.5V OPHTHALMOSCOPE
LOHR	3.5V RECHARGEABLE HANDLE ONLY
LODCD	DESK CHARGER DOUBLE CELL WITH TRANSFORMER
LOT	POWER ADAPTER / TRANSFORMER ONLY

#### 6 Standard Apertures



Cobalt blue filter, Micro spot, small spot, large spot, Slit, Fixation aperture.

#### With Red-Free Filter



Cobalt blue filter, Micro spot, small spot, large spot, Slit, Fixation aperture.



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107  
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134  
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au  
File Name: LOTPCSET\_Diagnostic-Set\_v4-0

First Prepared: 19/11/14 Eevien Tan  
Release Approval: 19/11/14 Lee Carter  
Version Control: V4.0 Page 1 | 1