

Mini Otoscope/Ophthalmoscope Pocket Diagnostic Set



The Liberty pocket or mini Otoscope/Ophthalmoscope diagnostic set is super convenient, portable and value-packed with strong metal pocket clips for convenient carrying while in use. The ophthalmoscope has 6 diagnostic beams and the closed head fibre optic otoscope comes complete with 4 extra speculae (ear tips). Use both in the consulting room and when on the move.

Mini Fibre Optic Otoscope features

- The Liberty Fibre Optic Mini Otoscope is a new and improved design, with a new switch and new modified bayonet fitting for a more secure attachment of speculae
- When the 2.5x magnifying glass is removed from the otoscope, external instruments can be safely inserted
- Otoscope has a bright light 2.5 volt F.O. light source and battery handle (uses 2 x AA size batteries which are included)

Mini Ophthalmoscope features

- Halogen light for true tissue colour and consistent, long-lasting illumination
- Apertures include large and small spot, semi-circle and red-free filter
- Rubber brow rest prevents scratching of eye-glasses
- Dust protection for a maintenance free life
- Pocket clip with integral on/off slide switch for ambidextrous use and battery handle (uses 2 x AA size batteries which are included)

Supplied in a solid protective carry case



Code	Description
LOTOPSET	LIBERTY MINI OTOSCOPE/OPHTHALMASCOPE DIAGNOSTIC SET
LOTETD2.5	DISPOSABLE EAR TIPS 2.5MM (PK 50)
LOTETD4	DISPOSABLE EAR TIPS 4MM (PK 50)

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: LOTOPMSET_Diagnostic-Set-Mini_v4-0



First Prepared: 19/11/14 Eevien Tan
Release Approval: 19/11/14 Lee Carter
Version Control: V4.0 Page 1 | 1