

Otoscope/Retinoscope Set

Liberty Rechargeable Diagnostics



The Liberty Fibre Optic Otoscope 3.5V has a sealed system for pneumatic otoscopy and accurate diagnosis of otitis media and a swivelling wide-angle viewing lens.

Features

- 3.5V Fibre Optic Otoscope
- Rechargeable handle, 120 minutes run time
- Swivelling wide-angle magnifying lens
- Reusable (specula) ear tips
- Compatible with Welch Allyn disposable ear tips

The Liberty Rechargeable Streak Retinoscope is designed exclusively to determine the refraction of the eye. The streak retinoscope is found by most practitioners to be easy to use, fast, accurate and especially valuable in determining the axis of astigmatism.

The Liberty Rechargeable Streak Retinoscope has multi-coated optics, a single control sleeve, totally dust proof, and a powerful bright bulb.

Features

- 3.5V Streak Retinoscope
- Rechargeable handle, 120 minutes run time
- Multi-coated optics
- A single control sleeve for vergence and rotation. Comfortable operation
- Totally dust proof, maintenance free
- 3.5V ABS handle-for long lasting operation
- Superior streak quality enables easier and quicker refractions-even in patients with small pupils
- The powerful, bright bulb ensures easy recognition of neutralisation
- The lowest position gives a divergent beam

WARRANTY:
 1 Year Manufacturers Warranty excluding bulbs and globes. The Rechargeable Handle (with 120 minutes run time) has a 2 Year Warranty. There is no internal battery. A new handle will need to be purchased once the handle no longer works.

Set Includes

- 1 x 3.5V Otoscope Head
- 1 x 3.5V Retinoscope Head
- 2 x Rechargeable handle, 120 minutes run time
- 1 set (2.75mm and 4.25mm, 5 pcs each) Reusable Specula (Ear Tips)
- 1 x Desk Charger Double Cell
- 1 x Power Adapter/Transformer



Order Code	Description
LOTRS	3.5V OTOSCOPE/RETINOSCOPE RECHARGEABLE SET
WA52432-U	DISPOSABLE EAR TIPS 2.75MM (PK 850)
WA52432-34	DISPOSABLE EAR TIPS 2.75MM (PK 34)
WA52434-U	DISPOSABLE EAR TIPS 4.25MM (PK 850)
WA52434-34	DISPOSABLE EAR TIPS 4.25MM (PK 34)
LOTCRSG	3.5V OTOSCOPE SPARE BULB
LOSRSG	3.5V STREAK RETINOSCOPE SPARE BULB
LOHR	3.5V RECHARGEABLE HANDLE ONLY
LODCD	DESK CHARGER DOUBLE CELL WITH TRANSFORMER
LOT	POWER ADAPTER/TRANSFORMER ONLY

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: LOTRS_Otoscope-Retinoscope_v4-0



First Prepared: 19/11/14 Eevien Tan
 Release Approval: 19/11/14 Lee Carter
 Version Control: V4.0 Page 1 | 1