

# Resuscitation Sets

## Adult, Child & Infant



Economical, durable silicone manual resuscitator designed to meet the need for high performance and low cost. Constructed from poly carbonate and silicone material that is autoclavable up to 134°C

### Adult Resuscitation Set

- Adult resuscitator 1500ml with 60cm H2O pop-off valve
- 2500ml Oxygen reservoir
- Two adult masks #5, #4
- Two guedel airways
- See through carry case with instruction manual



### Child Resuscitation Set

- Child resuscitator 550ml with 40cm H2O pop-off valve
- 2500ml Oxygen reservoir
- Two child masks #3, #2
- Two guedel airways
- See through carry case with instruction manual



### Infant Resuscitation Set

- Infant resuscitator 280ml with 40cm H2O pop-off valve
- 600ml Oxygen reservoir
- Two infant masks #1, #0
- Two guedel airways
- See through carry case with instruction manual



Code	Description
LRESSETA	ADULT MANUAL SILICONE RESUSCITATION SET
LRESSETC	CHILD MANUAL SILICONE RESUSCITATION SET
LRESSETI	INFANT MANUAL SILICONE RESUSCITATION SET

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107  
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134  
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au  
 File Name: LRESSET\_Resuscitation-Sets\_v2-0



First Prepared: 02/02/15 Eevien Tan  
 Release Approval: 02/02/15 Lee Carter  
 Version Control: V2.0 Page 1 | 1